

Protein *N*-myristoylation and prenylation in an insect cell-free protein synthesis system and their identification by mass spectrometry

Takashi Suzuki¹, Koko Moriya², Masaaki Ito¹, Toru Ezure¹, Masamitsu Shikata¹, Eiji Ando¹, Toshihiko Utsumi², Osamu Nishimura^{1,3}, Susumu Tsunasawa^{1,3}

1 Life Science Laboratory, Analytical and Measuring Instruments Division, Shimadzu Corporation, Kyoto, Japan
 2 Department of Biological Chemistry, Faculty of Agriculture, Yamaguchi University, Yamaguchi, Japan
 3 Institute for Protein Research, Osaka University, Osaka, Japan

Objective

The functional analysis of proteins in postgenomic studies has been attracting considerable attention, and interest in analyzing posttranslational modifications (PTMs) of proteins is increasing. Cell-free protein synthesis systems are assumed to be powerful tools for these objectives. Protein *N*-myristoylation and prenylation are the important lipid modifications of proteins, and they play crucial roles in regulating reversible protein-membrane and protein-protein interactions. Metabolic labeling is an effective strategy for analysis of these modifications. A rabbit reticulocyte lysate cell-free system has been widely utilized for the strategy, because it contains all the components involved in these protein modifications. However, metabolic labeling cannot be used to identify the exact location of the modification.

In this study, to establish a novel and effective strategy to analyze these protein modifications, human truncated gelsolin (tGelsolin) and rhoC, which are a natural *N*-myristoylated and geranylgeranylated protein respectively, were synthesized using an insect cell-free protein synthesis system. Synthesized proteins were analyzed by MALDI-TOF MS and MALDI-quadrupole ion trap (QIT)-TOF MS.

Materials and methods

Transdirect *insect cell*, which is based on the *Spodoptera frugiperda* 21 (SF21) extract, is a commercial product of Shimadzu (Kyoto, Japan). Epitope-tagged human tGelsolin and rhoC were synthesized using the insect cell-free protein synthesis system with or without the addition of a specific substrate for *N*-myristoylation or prenylation, such as myristoyl-CoA, farnesyl pyrophosphate, and geranylgeranyl pyrophosphate. The affinity purified proteins were reduced and alkylated, and then digested with trypsin or lysyl endopeptidase (Lys-C). The mass spectra and MS/MS spectra were acquired in reflectron positive ion mode with an AXIMA-CFR-plus MALDI-TOF MS instrument and an AXIMA-QIT MALDI-QIT-TOF hybrid mass spectrometer (Shimadzu/Kratos, Manchester, UK).

Results

In the case of the tGelsolin, when myristoyl-CoA was added to the translation reaction mixture, a peak corresponding to the *N*-myristoylated tryptic peptide was clearly observed, whereas this peak was hardly detected when myristoyl-CoA was not added. The peptide peak was identified as the *N*-myristoylated peptide fragment by tandem MS analysis using MALDI-QIT-TOF MS. We also constructed the Gly-2 to Ala (G2A) mutant, in which the *N*-myristoylation motif was disrupted, and analyzed by the same strategy. *N*-myristoylation did not occur on the G2A mutant, whereas this mutant was found to be *N*-acetylated after removal of the initiating Met.

In the case of the rhoC, it was found that it could serve as a substrate for both prenyltransferases in the presence of either farnesyl or geranylgeranyl pyrophosphate, whereas geranylgeranylation was only observed when both prenyl pyrophosphates were added to the translation reaction mixture. These results indicated that the rhoC protein was a better substrate for geranylgeranyltransferase I than for farnesyltransferase.


Conclusions

These results indicate that the insect cell-free protein synthesis system, as is the case with the rabbit reticulocyte lysate system, possesses *N*-myristoyltransferases and both prenyltransferases. In conclusion, a combination of the cell-free protein synthesis system with MS is an effective strategy to analyze protein *N*-myristoylation and prenylation.


References

- Suzuki, T., Ito, M., Ezure, T., Shikata, M., et al., N-Terminal protein modifications in an insect cell-free protein synthesis system and their identification by mass spectrometry. *Proteomics*, 6, 4486-4495 (2006)
- Suzuki, T., Ito, M., Ezure, T., Shikata, M., et al., Protein prenylation in an insect cell-free protein synthesis system and identification of products by mass spectrometry. *Proteomics*, 7, 1942-1950 (2007). This paper was selected as an "In This Issue" of the PROTEOMICS journal. The key points of the work have been available for download in the June edition of the PROTEOMICS podcast (<http://www.podcast.proteomics-journal.com>).

N-myristoylation




Expression and purification of tGelsolin




PMF analyses of the tryptic digests of tGelsolin proteins

N-terminal structures	Theoretical	Observed
Wild-type M G L G L S •• V E R	1846.9	ND
Wild-type G L G L S •• V E R	1715.9	1716.0
Wild-type M A L G L S •• V E R	1860.9	ND
Wild-type M A L G L S •• V E R	1926.1	1926.2
G2A mutant A L G L S •• V E R	1729.9	1729.8
G2A mutant A L G L S •• V E R	1940.1	ND


Fig. 3. Probable N-terminal structures of tryptic digests.




MS/MS analyses



Prenylation




Expression and purification of rhoC



PMF analyses of Lys-C digests of the rhoC proteins

C-terminal structures	Theoretical
R R R G C P I L	1027.6
R R R G C P I L	1174.8
R R R G C P I L	1242.8

Fig. 8. Probable C-terminal structures of Lys-C digests.



MS/MS analyses

