PharmaReport Program
Instruction Manual

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Read this manual carefully before using this product. Please keep this manual for future reference.

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
Instruction Manual
IRsolution
PharmaReport Program

Read this manual carefully before using this product.
Please keep this manual for future reference.

Analytical & Measuring Instruments Division

Copyright © Shimadzu Corporation 2003

All rights are reserved, including those to reproduce this publication or parts there of in any form without permission in writing from Shimadzu Corporation.

1. Information in this publication is subject to change without notice and does not represent a commitment on the part of vender.
2. Any errors or omissions which many have occurred in this publication despite the utmost care taken in its production will be corrected as soon as possible, but not necessarily immediately upon detection.
3. Note the Shimadzu does not have any obligation concerning the effects resulting from the application of the contents of this manual.

Please acknowledge being not able to answer for the influence of the result of operation by the content of this book.

For safety use

IRsolution PharmaReport Program is an optional program for IRsolution to create a report for identification test for Pharmacopoeia works. Please remember the following items when using this product.

1. Please DO NOT use this product for purposes other than what is was designed for.
2. For safety reasons, follow the procedure described in the instruction manual
3. Obey warnings or caution descriptions.
4. DO NOT modify the product by yourself.

This instruction manual defines the levels of warning as follows;

Warning	Death or serious injury may result is this warning is not followed.
Caution	Injury may result if this caution is not followed.
Notice	This notice with increase the understanding and usefulness of the software.

Guide for using the Instruction Manual

The manual of the PharmaReport program for the IRsolution software is one this book composition.

This instruction manual explains functions of the PharmaReport program function offered as an option function of the IRsolution software in this book. For general operations of the IRsolution, refer to the Instruction Manual for the IRsolution software.

Contents

Chapter 1	Outline of the PharmaReport Program	6
1-1	Outline	6
1-2	Environment to operate the PharmaReport Program	6
Chapter 2	Setup.....	7
2-1	Confirmation of number.....	7
2-2	Setup	7
2-3	Print template	8
2-4	Notes	8
2-4-1	Selection of spectrum	8
2-4-2	Printed Scan parameter.....	9
Chapter 3	How to use of the program.	10
3-1	Data which should be prepared	10
3-2	Beginning of program	10
3-3	Outline of the operation.....	11
3-4	Load Standard (standard spectrum) button	12
3-5	Load Sample (sample spectrum) button	12
3-6	Peak table file.....	13
3-6-1	Load Peak Table (reading of peak table file).....	13
3-6-2	Edit Peak Table (making and edit of peak table file)	13
3-7	Calculation (calculation execution) button.....	15
Chapter 4	How to use of the program.	18
4-1	Data which should be prepared	18
4-2	Beginning of program	18
4-3	Outline of the operation	19
4-4	Operation method.....	19
Chapter 5	How to use of the program.	22
-	Peak detection function -	22
5-1	Data which should be prepared	22
5-2	Beginning of program	22
5-3	Outline of the operation.....	23
5-4	Load Sample (sample spectrum) button	23
5-5	Peak table file.....	24
5-5-1	Load Peak Table (reading of peak table file).....	24
5-5-2	Edit Peak Table (making and edit of peak table file)	25
5-6	Detection (peak detection) button	26
Chapter 6	How to use of the program.	28
-	Alteration check program -	28
6-1	Operation method.....	28

Chapter 1 Outline of PharmaReport Program

1-1 Outline

The PharmaReport Program for IRsolution is supplied as an option program which operates on the IRsolution software. The following functions are included in this program.

1. Report Creating

Creates an analysis report to compare a sample spectrum with a standard spectrum based on identification test in the pharmacopoeias. Peak wavenumbers and peak ratios are calculated and the program judges Pass (OK)/Fail (NG) based on the errors calculated Peak wavenumbers and Peak ratios.

2. Spectrum Comparison

Creates a report suitable to compare standard and sample spectra based on identification test in the pharmacopoeias.

3. Peak Detection

Detect peaks being described in identification test in the pharmacopoeias, and creates a report.

1-2 Environment to operate the PharmaReport Program

1. Version of IRsolution

The PharmaReport Program operates by IRsolution Version 1.10 or later. Please contact to you Shimadzu representatives to upgrade your IRsolution version, when you use older version of IRsolution.

2. Macro platform

The PharmaReport Program operates by using the Macro Platform of the IRsolution software. If the Macro Platform is not installed, the PharmaReport does not run. Please purchase and install the Macro Platform when it is not possession.

Chapter 2 Setup

2-1 Confirmation of parts

Check up necessary parts are included in your package before installation.

Name of articles	P/N	No.
Installation disk for PharmaReport	206-72488-91	1 set
Instruction Manual (this book)	206-95130	1

2-2 Setup

This paragraph describes the procedure to install the PharmaReport Program to the IRsolution software.

Please confirm the IRsolution software and Macro Platform has been correctly installed in the computer before the installation.

1. Insert the installation disk in the floppy disk drive.
2. Select "Run" on the start button menu of Windows.
3. Input to the file name as "A:setup", and then clicks < OK > button.
4. The installation program starts.
5. Click < OK > button, and the installation is finished.
6. "ParmReport.lcm" is copied to "... ¥Program files¥Shimadzu¥IRsolution¥macros" folder. "Kyokuhou3.ptm", "Kyokuhou2.ptm", "Kyokuhou1NG.ptm", and "Kyokuhou1OK.ptm" are copied onto "... ¥Program files ¥Shimadzu ¥IRsolution ¥PrintTemplates" folder.

2-3 Print template

When IRsolution is not installed in “C:\Program files\Shimadzu\IRsolution”, it is necessary to correct the print template file. Please correct four files of “Kyokuhou3.ptm”, “Kyokuhou2.ptm”, “Kyokuhou1NG.ptm”, and “Kyokuhou1OK.ptm” as follows.

1. Run the “Report” function, and click < No > button on the print screen once. As a result, “...\Program files\Shimadzu\IRsolution\data\kyokuhou\report.txt” and “report2.txt” are created.
2. Select the [Print] tab on IRsolution.
3. Select [File]-[Open] and open “Kyokuhou3.ptm”.
4. Select the evaluation result text frame.
5. Specify the “Report.txt” or “Report2.txt” to the result text frame on the “RTF Filename” column in object property window.
6. Select [File]-[Save].
7. Other “Kyokuhou2.ptm”, “Kyokuhou1NG.ptm”, and “Kyokuhou1OK.ptm” are similarly corrected.

2-4 Notes

Please note the following points when you use this program.

2-4-1 Selection of spectrum

In this program, data with an “active pointer” in the selected spectrum file is used to calculation.

In IRsolution, when various Data manipulation such as Baseline correction, Kubelka-Munk conversion, etc. are done to the measured data, the Data manipulation result is saved as a “Child” data of the original data in the has the tree-structured file. Data with an “Active pointer” is displayed first when a file consisted of two or more data is opened. In this program, data with an “active pointer” in the selected spectrum file is used to calculation. An active pointer is usually set to the data treated at the end. When an active pointer is not set to target data, the file is opened with IRsolution, and mouse’s right button is clicked on the tree view. Then select an active pointer of the displayed context menu, and close the file.


Figure 2-1 Tree view and Active pointer

2-4-2 Scan parameters on Report

When the Raw (scanned) data is used by this program, Scan parameters are printed on the report. However, when an active pointer is set to manipulated data, the Scan parameters are not printed, and leave blank. The reason is that Manipulated data does not include scan parameters because it is not raw data.

Chapter 3 How to use of the program. - Report creating function -

The Report function creates an analysis report to compare a sample spectrum with a standard spectrum based on identification test in the pharmacopoeias. Peak wavenumbers and peak ratios are calculated and the program judges Pass (OK)/Fail (NG) based on the errors calculated Peak wavenumbers and Peak ratios.

3-1 Data which should be prepared

- Standard spectrum (Transmittance)
- Sample spectrum (Transmittance)
- Peak table (It is made in this program.)

Note

- ◆ In this program, data with an “active pointer” in the selected spectrum file is used to calculation.
- ◆ When the Raw (scanned) data is used by this program, Scan parameters are printed on the report. However, when an active pointer is set to manipulated data, the Scan parameters are not printed, and leave blank.

3-2 Starting of program

1. Start IRsolution.
2. Click the [Macro] tab to display the macro platform screen.
3. Select [Macro]-[Open Macro Container] to open “...¥Program Files ¥Shimadzu ¥IRsolution¥Macros¥PharmaReport.lcm”. The program is displayed in the macro toolbar.
4. Select “PharmaReport” under the “PharmaReport.lcm”, and then select [Macro]-[Run] to execute the program.

Note The program does not start even if the  button in the macro toolbar is clicked.

5. The dialog box in Figure 3-1 opens.


Figure 3-1. Selection of program

Button	Content
Printout	Runs the Spectrum Comparison function.
Peak Detection	Runs the Peak detection function.
Report	Runs the Report creating function.

6. Then, click < Report > to open the dialog box in Figure 3-2, and the Report creating function for Pharmacopoeias starts.


Figure 3-2. "Report creating function" screen

Button	Content
Load Standard	Loads a standard spectrum.
Load Sample	Loads a sample spectrum.
Load Peak Table	Loads a Peak table file which has already been made.
Edit Peak Table	Creates/Edits a Peak table file.
Calculation	Calculates and creates a Report.
Exit	Closes the program.

3-3 Outline of the operation

The outline of the operation of this program is as follows.

1. Load a standard spectrum with < Load Standard > button.
2. Load a sample spectrum with < Load Sample > button.
3. Load a Peak table file which has already been made with < Load Peak Table > button, or create a Peak table file with < Edit Peak Table > button.
4. Click the < Calculation > button to create a Report.

3-4 <Load Standard> Spectrum button

It is a button to prepare a standard spectrum file necessary for the calculation.

1. Click the < Load Standard > button to open the “Load Standard Spectrum” dialog box.
2. Select the standard spectrum
3. Clicking the <Open > button returns to the dialog box in Figure 3-2.


Figure 3-3. "Load Standard Spectrum" dialog box

3-5 <Load Sample> Spectrum button

It is a button to prepare a sample file necessary for the calculation.

1. Click the < Load Sample > button to open the “Load Sample Spectrum” dialog box.
2. Select the sample spectrum.
3. Clicking the <Open > button returns to the dialog box in Figure 3-2.


Figure 3-4. "Load Sample Spectrum" dialog box

3-6 Peak table file

Up to 10 peaks can be specified to detect. Up to 4 peaks can be specified to calculate the Transmittance Ratio to evaluate the intensity. The set of the peak wavenumber is saved as a data file.

The peak table file can be created by clicking < Load Peak Table > button in Figure 3-2, or edited by clicking < Edit Peak Table > button.

3-6-1 <Load Peak Table> button

The peak table file which has already been made can be loaded by clicking <Load Peak Table> button in Figure 3-2.

1. The dialog box in Figure 3-5 opens.


Figure 3-5. Load Peak File dialog box

2. The extension of the peak table file is “.wdtr”.
3. Select a Peak table file and then, click the <Open> button to return to the dialog box in Figure 3-2.
4. Clicking the < Edit Peak Table > button explains next is used to create the Peak table newly.

3-6-2 <Edit Peak Table> button

The peak table file can be created or edited by clicking < Edit Peak Table > button in Figure 3-2.

1. The dialog box in Figure 3-6 opens.

Use	No.	Wave No.	Tolerance(cm-1)
<input checked="" type="checkbox"/>	1	540	5
<input checked="" type="checkbox"/>	2	906	5
<input checked="" type="checkbox"/>	3	1028	5
<input checked="" type="checkbox"/>	4	1155	5
<input checked="" type="checkbox"/>	5	1370	5
<input checked="" type="checkbox"/>	6	1585	5
<input checked="" type="checkbox"/>	7	1600	5
<input checked="" type="checkbox"/>	8	1942	5
<input checked="" type="checkbox"/>	9	2850	5
<input checked="" type="checkbox"/>	10	3060	5

Evaluation Peak No. 4

A 1 B 3 C 7 D 10

Tolerance(0-1) 1+/- 0.3

Figure 3-6. "Edit Peak Information" dialog box

2. Input peak wavenumbers into "Wave No." column, and the range of the Tolerance of the peak wavenumber into the "Tolerance" column. It is possible to input up to ten peaks. The check mark is put to "Use" columns of the wavenumber used to calculate. Even if the wavenumber is input, the one without check mark is not used to calculate.
3. First of all, number of peaks to be used for Peak ratio evaluation to "Evaluation Peak No.". Then input the numbers of the peak wavenumber to be used for evaluation. It is possible to specify up to four peaks. The range of the Tolerance compared with the peak is input to the "Tolerance" column.
4. When the input ends, < OK > button is clicked. The dialog box in Figure 3-7 is displayed, and the created peak table file is saved. The extension of the peak table file is ".wdtr".
5. Clicking the <Save> button returns to the dialog box in Figure 3-2.


Figure 3-7. "Save Peak File" dialog box

3-7 <Calculation> button

The calculation is executed, and judged. The result is displayed on the screen and printed.

1. The <Calculation> button will be activated when all of Loading Peak Table, Loading Sample file and Loading Standard file are done.
2. When < Calculation > button is clicked, the calculation is executed according to the following procedures.
 - (1) The peak of standard spectrum is detected referring to the peak table.
 - (2) The peak of the sample spectrum is detected referring to the detected peak list of standard spectrum.
 - (3) Errors (differences of peak wavenumber of standard from one pf sample) are calculated.
 - (4) Peak ratios in Transmittance are calculated for the standard spectrum.
 - (5) Peak ratios in Transmittance are calculated for the sample spectrum.
 - (6) Ratio of peak ratios in (5) to peak ratios in (4) for evaluation.
 - (7) Prints out the result of (3) and (6). At this time, Judgment becomes "NG" when one of errors of peak wavenumber or one of errors of peak ratio is not satisfied Tolerance.
 - (8) When all of all listed errors are satisfied Tolerances, Judgment is "OK".
3. The spectrum and the Judgment ("OK" "NG") are displayed.

Figure 3-8. Evaluation

4. The dialog box in Figure 3-9 opens when < OK > button is clicked.

Figure 3-9. Print dialog box

5. Click <Yes> button to printout a Report such as Figure 3-10. The program ends if <No> button is clicked.


Figure 3-10. Printed report example

The report consists of the Scan parameter, Judgment (OK or NG), spectra of standard and sample, and calculated table.

Markers to show the selected peaks for detection and evaluation are displayed.

The upper part of the calculated table shows the detected peak wavenumbers and their transmittance for both sample and standard spectra. The right row shows error which is calculated by “[Peak wavenumber of sample] – [Peak wavenumber of standard]”.

The lower part shows Peak ratios of evaluation peaks for both sample and standard spectra. The right row shows a error which is calculated by “[Peak ratio of sample] / [Peak ratio of standard]”.

Chapter 4 How to use of the program. - Spectrum comparison function -

The Spectrum comparison function creates a report suitable to compare standard and sample spectra based on identification test in the pharmacopoeias.

4-1 Data which should be prepared

- Standard spectrum (Transmittance)
- Sample spectrum (Transmittance)

Note

- ◆ In this program, data with an “active pointer” in the selected spectrum file is used to calculation.
- ◆ When the Raw (scanned) data is used by this program, Scan parameters are printed on the report. However, when an active pointer is set to manipulated data, the Scan parameters are not printed, and leave blank.

4-2 Starting of program

1. Start IRsolution.
2. Click the [Macro] tab to display the macro platform screen.
3. Select [Macro]-[Open Macro Container] to open “...¥Program Files ¥Shimadzu ¥IRsolution¥Macros¥PharmaReport.lcm”. The program is displayed in the macro toolbar.
4. Select “PharmaReport” under the “PharmaReport.lcm”, and then select [Macro]-[Run] to execute the program.

Note The program does not start even if the  button in the macro toolbar is clicked.

5. The dialog box in Figure 4-1 opens.


Figure 4-1. Selection of program

Button	Content
Printout	Runs the Spectrum Comparison function.
Peak Detection	Runs the Peak detection function.
Report	Runs the Report creating function.

- Then, click < Printout > button.

4-3 Outline of the operation

The outline of the operation of this program is as follows.

- Load a standard spectrum.
- Load a sample spectrum.
- The report is printed.

4-4 Operation method

- Click < Printout > button in Figure 4-1.
- The dialog box in Figure 4-2 opens.


Figure 4-2. "Load Standard Spectrum" dialog box

- Select a standard spectrum, and then click <Open > button.
- The dialog box in Figure 4-3 opens.


Figure 4-3. "Load Sample Spectrum" dialog box

5. Select a sample spectrum, and then click <Open > button.
6. The dialog box in Figure 4-4 opens.


Figure 4-4. Print dialog box

7. Click <Yes> button to printout a Report such as Figure 4-5. The program ends if <No> button is clicked.


Figure 4-5. Printed report example

Chapter 5 How to use of the program. - Peak detection function -

The Peak detection function detect peaks being described in identification test in the pharmacopoeias, and creates a report.

5-1 Data which should be prepared

- Sample spectrum (Transmittance)
- Peak table (It is made in this program.)

Note

- ◆ In this program, data with an “active pointer” in the selected spectrum file is used to calculation.
- ◆ When the Raw (scanned) data is used by this program, Scan parameters are printed on the report. However, when an active pointer is set to manipulated data, the Scan parameters are not printed, and leave blank.

5-2 Starting of program

1. Start IRsolution.
2. Click the [Macro] tab to display the macro platform screen.
3. Select [Macro]-[Open Macro Container] to open “...¥Program Files ¥Shimadzu ¥IRsolution¥Macros¥PharmaReport.lcm”. The program is displayed in the macro toolbar.
4. Select “PharmaReport” under the “PharmaReport.lcm”, and then select [Macro]-[Run] to execute the program.

Note

The program does not start even if the  button in the macro toolbar is clicked.

5. The dialog box in Figure 5-1 opens.


Figure 5-1. Selection of program

Button	Content
Printout	Runs the Spectrum Comparison function.
Peak Detection	Runs the Peak detection function.
Report	Runs the Report creating function.

6. Then, click < Peak Detection > to open the dialog box in Figure 5-2, and the Peak detection function for Pharmacopoeia starts.


Figure 5-2. "Peak detection function" screen

Button	Content
Load Sample	Loads a sample spectrum.
Load Peak Table	Loads a Peak table file which has already been made.
Edit Peak Table	Creates/Edits a Peak table file.
Detection	Detects peak and creates a Report.
Exit	Closes the program.

5-3 Outline of the operation

The outline of the operation of this program is as follows.

1. Load a sample spectrum with < Load Sample > button.
2. Load a Peak table file which has already been made with < Load Peak Table > button, or create a Peak table file with < Edit Peak Table > button.
3. Click the < Detection > button to create a Report.

5-4 <Load Sample> Spectrum button

It is a button to prepare a sample file necessary for the calculation.

1. Click the < Load Sample > button to open the "Load Sample Spectrum" dialog box.

2. Select the sample spectrum.
3. Clicking the <Open > button returns to the dialog box in Figure 5-2.


Figure 5-3. "Load Sample Spectrum" dialog box

5-5 Peak table file

Up to 10 peaks can be specified to detect. Up to 4 positions can be specified to remark. The set of the peak wavenumber is saved as a data file.

The peak table file can be created by clicking < Load Peak Table > button in Figure 5-2, or edited by clicking < Edit Peak Table > button

5-5-1 <Load Peak Table> button

The peak table file which has already been made can be read by clicking <Load Peak Table> button in Figure 5-2.

1. The dialog box in Figure 5-4 opens.


Figure 5-4. "Load Peak Table" dialog box

2. The extension of the peak table file is “.wdtp”.
3. Select a Peak table file and then, click the <Open> button to return to the dialog box in Figure 5-2.
4. Clicking the < Edit Peak Table > button explains next is used to create the Peak table newly.

5-5-2 <Edit Peak Table> button

The peak table file can be created or edited by clicking < Edit Peak Table > button in Figure 5-2.

1. The dialog box in Figure 5-5 opens.


Figure 5-5. “Edit Peak W. No.” dialog box

2. Input peak wavenumbers into “Wave No.” column. It is possible to input up to ten peaks. The check mark is put to “Use” columns of the wavenumber used to calculate. Even if the wavenumber is input, the one without check mark is not used to calculate.
3. First of all, number of positions to be used for Remark. Then input the wavenumbers for Remarks. It is possible to specify up to four positions.
4. When the input ends, < OK > button is clicked. The dialog box in Figure 5-6 is displayed, and the created peak table file is saved. The extension of the peak table file is “.wdtp”.
5. Clicking the <Save> button returns to the dialog box in Figure 5-2.


Figure 5-6. "Save Peak File" dialog box

5-6 <Detection> button

The peak detection is executed. The result is printed.

1. The <Detection> button will be activated when all of Loading Peak Table and Loading Sample file are done.
2. The dialog box in Figure 5-7 opens.


Figure 5-7 print dialog box

3. Click <Yes> button to printout a Report such as Figure 5-8. The program ends if <No> button is clicked.

The procedure from which the calculation is executed is as follows.

- (1) The peak of the sample spectrum is detected referring to the peak table.
- (2) When the peak is not detected, the error message is actually displayed.
- (3) The peak position is plotted on the screen.
- (4) When <detection> is selected, a spectrum, scan parameters, and peak information are printed such as Figure 5-8.
- (5) Close the program clicking "Cancel" button once, and select [File]-[SaveAs] command to save the spectrum with peak marks.

Figure 5-8. Printed report example

Chapter 6 How to use of the program. - Alteration check program -

The Alteration check function is included in this program. The Alteration check function is a program which checks whether this program is illegally altered.

6-1 Operation method

1. Select Alteration check from [Start]-[Program].
2. When “Shimadzu user authentication tool” is installed, Login dialog box of Figure 6-1 is displayed. Input User ID and password. When the tool is not installed, the dialog in Figure 6-2 is displayed.

Figure 6-1. “Log in” screen

3. “Alteration check program” screen in Figure 6-2 is displayed.


Figure 6-2 “Alteration check program” screen

4. Select “PharmaReport program”, and click < check > button.
5. When the screen in Figure 6-3 is displayed, click the < execution > button.


Figure 6-3. "Alteration check" screen

6. The Alteration check result is displayed in the dialog.


Figure 6-4. Check result

7. When < Note pad > button is clicked, the check result is displayed. It is possible to print from this screen. The check result cannot be corrected.


Figure 6-5. Check result

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Copyright © 2002 株式会社 島津製作所