

Simon Ashton<sup>1</sup>; Neil Loftus<sup>1</sup>; Chris Titman<sup>1</sup>; Albert Koulman<sup>2</sup>  
<sup>1</sup>Shimadzu, Manchester , United Kingdom ; <sup>2</sup>Lipidomics Biomarker Research Group, Elsie Widdowson laboratory, MRC, Cambridge, UK

Introduction

Untargeted metabolite profiling has become an integral part of systems biology research and is helping to further our understanding of the physiological and patho-physiological processes related to diet, environmental factors, disease and the impact of pharmaceutical treatments. In practical terms global metabolite profiling needs to take into account not only the inherent complexity and diversity of human metabolism but also analytical technologies that deliver robust, reproducible and high quality data. To help standardize and objectively measure metabolite levels in large sample numbers a standard biological Quality Control has become a widely accepted tool in metabolic profiling to assess system performance and to filter ion signals which show idiosyncratic response throughout the batch analysis. This paper discusses the application of an ion trap-time of flight mass spectrometer to measure endogenous metabolite levels in human plasma profiling studies.

Methods

The plasma of healthy volunteers was used, either undiluted, with 20% Ringer addition, with glucose in Ringer addition, with palmitate in Ringer addition, with N-octanulphosphingone in Ringer or with glucose, palmitate and N-octanulphosphingone together in palmitate. This resulted in sample set of 60 different samples. By combining aliquots of all samples a pooled QC samples was created. Plasma samples were analysed by Prominence HPLC coupled with an electrospray ion trap-time of flight mass spectrometry (LCMS-IT-TOF, Shimadzu Corporation, Kyoto, Japan) using a 0.1% formic acid / acetonitrile gradient with a sample cycle time of less than 10 minutes.




Figure 1. Mass chromatograms for a series of phospholipids from a pooled QC human plasma sample. Typical chromatographic peak width: 5-7 seconds; MS scan speed: 100msc; m/z 170-900; external mass calibration; ion accumulation time: 20msec. In MS mode 10 spectra per second were acquired. In polarity switching mode, the switching time was 100msec.

Results

In metabolite profiling studies it is important to assess the system reproducibility throughout the batch analysis. In this study the pooled QC sample was used to characterize the reproducibility of the 'system'. Whilst the FDA suggests that variability of  $\pm 15\%$  of the nominal value represents an acceptable degree of reproducibility, in long term profiling studies the tolerance is often relaxed to between 20-40%. In this study, ion signals which resulted in a relative standard deviation (RSD) <20% and retention time RSD<1% were considered in subsequent principal component analysis.

Lipid ID	Theoretical m/z (M+H) <sup>+</sup>	Avg RT Mins	Min RT	Max RT	%RSD RT	%RSD Area	Ion signal intensity [Average response]
GPCho(14:0/0:0)	468.3085	2.395	2.392	2.403	0.13	2.26	2,173,959
GPCho(18:1/0:0)	480.3085	2.650	2.654	2.665	0.13	4.38	685,760
GPCho(0:16:1)	480.3449	2.654	2.655	2.666	0.13	5.54	780,277
GPCho(18:0/0:0)	482.3241	2.490	2.482	2.493	0.13	3.06	1,321,240
GPCho(0:16:0)	482.3605	2.490	2.482	2.493	0.13	3.81	649,604
GPCho(16:1/0:0)	494.3241	2.451	2.451	2.462	0.13	2.33	4,722,352
GPCho(16:0/0:0)	496.3398	2.541	2.579	2.590	0.13	0.93	50,235,899
GPCho(20:0/0:0)	502.2928	2.520	2.518	2.529	0.12	4.61	594,206
GPCho(0:18:1)	508.3762	2.639	2.707	2.718	0.12	4.83	423,688
GPCho(18:1/0:0)	518.3241	2.405	2.407	2.417	0.13	3.62	1,241,123
GPCho(18:2/0:0)	520.3398	2.509	2.507	2.518	0.13	0.76	36,334,718
GPCho(18:1/0:0)	522.3554	2.640	2.634	2.645	0.12	1.48	25,727,526
GPCho(18:0/0:0)	524.3711	2.780	2.778	2.789	0.11	1.91	18,414,007
GPCho(22:5/0:0)	526.2938	2.654	2.501	2.512	0.13	9.65	272,533
GPCho(20:5/0:0)	542.3241	2.507	2.507	2.519	0.14	8.75	1,206,508
GPCho(20:0/0:0)	544.3398	2.512	2.507	2.518	0.13	1.02	6,262,232
GPCho(26:0/0:0)	546.3554	2.575	2.573	2.585	0.13	2.72	1,846,463
GPCho(22:5/0:0)	568.3398	2.492	2.486	2.497	0.13	2.35	2,762,070
SM(d18:1/24:0)	675.5436	5.550	5.578	5.592	0.10	2.23	1,838,330
SM(d18:1/15:0)	689.5592	3.870	3.692	3.702	0.10	8.05	132,281
SM(d18:1/16:1)	701.5592	3.618	3.624	3.638	0.10	2.71	1,477,189
GPCho(0:34:3)	742.5745	4.030	4.043	4.057	0.09	2.23	1,307,432
GPCho(0:34:2)	744.5902	4.066	4.054	4.069	0.09	4.59	946,724
GPCho(34:4)	754.5391	3.699	3.714	3.729	0.09	4.07	548,470
GPCho(34:3)	755.5538	3.758	3.767	3.781	0.10	2.37	5,074,504
GPCho(34:2)	758.5694	3.937	3.931	3.946	0.09	0.64	43,017,246
GPCho(34:1)	760.5851	4.095	4.096	4.111	0.09	1.05	19,206,390
GPCho(36:6)	764.5589	3.877	3.876	3.891	0.10	3.27	473,322
GPCho(36:5)	765.5745	4.017	4.022	4.037	0.09	2.58	279,456
GPCho(36:4)	767.5902	3.927	4.092	4.106	0.09	5.52	572,722
GPCho(36:3)	780.5538	3.768	3.764	3.779	0.10	2.12	15,194,035
GPCho(36:2)	782.5694	3.907	3.911	3.926	0.10	1.22	12,274,519
GPCho(36:1)	784.5851	3.977	3.980	3.995	0.09	1.31	29,450,182
GPCho(38:2)	786.6007	4.145	4.149	4.164	0.09	1.71	20,403,293
GPCho(0:38:7)	790.5745	3.937	3.953	3.968	0.10	7.66	30,440
GPCho(38:6)	792.5902	4.533	4.539	4.555	0.09	5.89	1,332,860
GPCho(38:5)	794.6058	4.075	4.074	4.088	0.09	3.13	2,951,349
GPCho(38:4)	796.6215	4.265	4.259	4.275	0.09	7.74	636,674
GPCho(38:3)	804.5538	4.007	4.011	4.028	0.10	6.76	420,325
GPCho(38:2)	806.5694	3.847	3.850	3.865	0.10	1.72	26,501,408
GPCho(38:1)	808.5851	3.996	3.974	3.990	0.10	3.07	7,476,945
GPCho(38:0)	810.6007	4.125	4.134	4.149	0.09	2.69	6,308,132
GPCho(38:3)	812.6164	4.205	4.206	4.221	0.09	5.03	4,405,377
GPCho(0:40:6)	820.6215	4.195	4.087	4.090	0.14	8.46	210,117
GPCho(40:7)	832.5851	3.897	3.889	3.905	0.10	3.75	1,548,653
GPCho(40:6)	834.6007	4.065	4.064	4.080	0.10	2.87	5,463,972
GPCho(40:5)	836.6164	4.125	4.066	4.140	0.08	4.34	1,374,137

Table 1. The table above shows the reproducibility of the system to a number of phospholipids with different signal intensities for a pooled QC sample (n=17; 1 $\mu$ L injection volume; interspersed between every 12 samples; total run sequence time was approximately 25 hours; [batch analysis size was 129 files which included 17 QC samples]). In this analysis the ion signal reproducibility was less than 10% for all phospholipids measured and the retention time variation was equal to or less than 0.1%. Peak areas calculated using LabSolutions software.

Results




Figure 2. Principal Component Analysis. Pooled QC data are closely clustered together and show no 'run order' change in signal response over the analytical run. [The data presented above consider the 1 $\mu$ L injection volume data only].

Measured	Theoretical	Mass accuracy	RT Range	Ion RT	RT #	#	Aligned data array - filtered - sorted according to ion intensity	q07
520.3398	520.3398	1.73	002.44 - 002.53	2.507	17	112	25833396 27191868 25643340 26505175 23366279 26324919 26001039	q05
408.3857	408.3834	5.63	003.18 - 003.25	3.231	17	112	24519187 24294177 24475130 24529375 23541062 23716732 24353098	q06
758.5709	758.5694	1.98	003.90 - 003.96	3.933	17	112	22881968 22515941 23098426 22344578 22928781 22380036 22200799	q07
522.3548	522.3554	-1.15	002.56 - 002.65	2.654	17	112	20680254 19967149 19829518 20397499 19382161 19912099 19780754	q08
524.3598	524.3711	-2.48	002.77 - 002.79	2.777	17	112	15444494 16208490 15845183 16223458 15951256 16068710 16300140	q09
820.5467	820.5538	-0.37	003.71 - 003.87	3.854	17	112	15449731 15082963 14747676 15384993 15073221 14458811 14470462	q10
786.6007	786.6007	-0.25	004.09 - 004.17	4.151	17	112	12328697 13171423 12790519 13326723 12902928 13523591 12646797	q11
782.5695	782.5694	0.13	003.78 - 003.83	3.914	17	112	12207563 11945942 12216128 11906726 11927999 11796744 12129750	q12
784.5853	784.5851	0.25	003.95 - 004.01	3.982	17	112	9879755 9910688 10410558 10107200 9995321 10042110 9924339	q13
780.5538	780.5528	1.28	003.68 - 003.79	3.767	17	112	9796279 9787057 10139088 9957328 9580785 9660253 9628606	q14

Table 2. Following the alignment of the raw data files the data array was filtered to include ions above a set tolerance (%RSD area <20%; %RSD RT <2%; minimum number of pooled QC ions was >80%). Plasma samples prepared by protein precipitation with cold acetonitrile typically results in a data array dominated by phospholipid ion signals. The table shows the average mass accuracy measured throughout the batch analysis (the mass accuracy was calculated for each sample [8s] using the m/z value for the most intense ion in the spectrum bin). For example, the average mass accuracy throughout the run of 112 samples and 17 QC controls for m/z 520.3398 was 1.73ppm (n=129).

Conclusions

- This study reports the application of a LCMS-IT-TOF to human plasma metabolite profiling. Using a scan speed of 10 scans per second the peak area variance for a number of phospholipid ion signals was less than 10% with the retention time variance less than or equal to 0.10% for peak widths typically between 5-7 seconds.
- Following the alignment of ion signals and applying filtering parameters (ion signals present in the pooled QC must vary less than 20% in peak area throughout the run and must be present in 80% of the QC samples) 916 ions were submitted for statistical analysis using Simca-P (Umetrics) . The PCA plot shows a cluster of QC signals indicating robust detection of the ion signals detected.