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糖	アミノ酸	その他 (有機酸、アミンなど)
Arabitol	2-Aminobutyric acid	2-Aminoethanol
Erythritol	2-Aminopimelic acid	2-Hydroxyglutaric acid
Erythrulose	3-Aminoglutaric acid	2-Ketobutyric acid
Fructose	3-Aminoisobutyric acid	2-Ketoglutaric acid
Fructose	3-Aminopropanoic acid	3-Hydroxybutyric acid
Galactose	4-Aminobutyric acid (GABA)	3-Hydroxyisovaleric acid
Galacturonic acid	5-Oxoproline	3-Hydroxypropionic acid
Glucaric acid	5-Aminovaleric acid	Aconitic acid
Glucosamine	Alanine	Allantoin
Glucose	Anthranilic acid	Benzoic acid
Glucuronic acid	Arginine	Caproic acid
Glyceric acid	Asparagine	Citramalic acid
Glycerol	Aspartic acid	Citric acid
Inositol	Glutamic acid	Fumaric acid
Lyxose	Glycine	Glycolic acid
Maltose-meto	Histidine	Hypotaurine
Psicose	Homoserine	Isocitric acid
Rhamnose	Isoleucine	Lactic acid
Ribulose	Leucine	Lauric acid
Trehalose	Lysine	Maleic acid
Xylose	Methionine	Malic acid
Xylulose	O-Acetylserine	Malonic acid
	Ornithine	Mesaconic acid
糖リン酸	Phenylalanine	Methylsuccinic acid
3-Phosphoglyceric acid	Proline	Niacinamide
Dihydroxyacetone phosphate	Saccharopine	Nicotinic acid
Fructose 6-phosphate	Serine	O-Phosphoethanolamine
Glucose 6-phosphate	Threonine	Oxalic acid
Glycerol 3-phosphate	Tyrosine	Palmitic acid
Mannose 6-phosphate	Tryptophan	Phosphoric acid
Ribose 5-phosphate	Valine	Protocatechuic acid Putrescine
Sedoheptulose 7-phosphate		Pyrogallol
	核酸	Pyruvic acid
	Adenine	Ribonolactone
	Adenosine	Stearic acid
	Adenosine monophosphate	Succinic acid
	Guanine	Tartaric acid
	Guanosine	Threonic acid
	Hypoxanthine	Tryptamine
	Inosine	Urea
	Uracil	Uric acid
	Uridine	
	Xanthosine	

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© Shimadzu Corporation, 2020