# Application News

## No. L529B

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- Mercaptopropionic Acid Reagent  
Add 10 μL of 3-mercaptopropionic acid into 10 mL of 0.1 mol/L borate buffer.
- OPA Reagent  
Add 0.3 mL of ethanol into 10 mg of *o*-phthalaldehyde and dissolve completely. Then add 0.7 mL of 0.1 mol/L borate buffer and 4 mL of pure water.
- Mercaptopropionic Acid / OPA solution  
Mix 300 μL of Mercaptopropionic Acid Reagent and 600 μL OPA Reagent.
- FMOC Reagent  
Add 10 mg of 9-fluorenylmethyl chloroformate into 50 mL of acetonitrile.
- Mobile phase A : 20 mmol/L Sodium acetate buffer (pH6)  
Add 2.67 g of sodium acetate trihydrate and 41 μL of acetic acid into 1000 mL of pure water.
- Mobile phase B : Water/Acetonitrile = 100/900
- Mobile phase C : 20 mmol/mL Sodium acetate buffer (pH5) containing 0.5 mmol/L EDTA-2Na  
Add 0.19 g of EDTA-2Na, 2.03 g of sodium acetate trihydrate and 308 μL of acetic acid into 1000 mL of pure water.
- phosphoric acid aqueous solution  
Add 0.5 mL of phosphoric acid into 100 mL of pure water.

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
1, Aspartic Acid 2, Glutamic Acid 3, Asparagine 4, Serine 5, Glutamine 6, Histidine 7, Glycine 8, Threonine 9, Arginine 10, Alanine 11, Tyrosine 12, Methionine 13, Valine 14, Cystine 15, Tryptophan 16, Phenylalanine 17, Isoleucine 18, Leucine 19, Proline 20, Lysine

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