Application News

No. B92

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View experiment results

Select a test from the list below (or the target on the left) to display the results of the test in greater detail.

Spot ID	Sample ID	Expected Mass	Actual Mass	Test Result
E2	Cal		Not found	Pass
E3	Exendin-4	4187.6	4187.55	Pass
E4	Exendin-4	4187.6	4187.85	Pass
D2	Exendin-4	4187.6	Not found	Fail
D3	Exendin-4	4187.6	4187.25	Query

QC Criteria

Select the criteria that constitutes a QC test. A QC test is marked as Pass if none of the Fail or Query criteria are met

Include in test	Criteria	Value	Measure	Fail	Query	Special Condit..
<input type="checkbox"/>	Mass tolerance is greater than	1	Da	<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	S/N threshold is less than	50		<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input type="checkbox"/>	Minimum resolution	1000		<input checked="" type="radio"/>	<input type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/>	Sum of adducts is greater than	50	%Area	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>
<input checked="" type="checkbox"/>	Maximum impurity level	50	%Area	<input type="radio"/>	<input checked="" type="radio"/>	<input type="radio"/>
Range Criteria						
<input type="checkbox"/>	S/N is within range			10	45	
<input type="checkbox"/>	Sum of adducts is within range		Apex(mV)	1	7	
<input type="checkbox"/>	Min > max acceptable impurity level		%Area	60	90	

Special Conditions - Optional additional specific criteria that may be applied

Mass NOT within tolerance but Adducts found: Fail Query

Define the Mass calculation Type: Monoisotopic Average


Define Adducts... 1 selected
 Define Amino Acids... 0 selected
 Define Matrix Ions... 0 selected
 Define Impurities... 0 selected

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b)

Detected Ions			
Ion Type	Description	Non-Ionised	Ionised
Target Adduct	OH	4202.27	4203.27
Impurities detected			
Non-Ionised	Ionised	Percent of target by area	Percent of target by intensity
4216.80	4217.80	146.88%	105.38%
2100.95	2101.95	10.78%	34.80%
4246.25	4247.25	5.39%	24.45%
Test criteria			
Criterion	Required value	Actual value	Result
Sum of adducts greater than/within range	50	557.45	Query
Maximum impurity level	50	163.04	Query

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