Application
Data Sheet

No. 103

GC-MS

Gas Chromatograph Mass Spectrometer

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Analysis of Metabolites Extracted from Human Embryonic Stem Cells
using GC-MS

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1	Adenine-2TMS	43	Glycine-2TMS	85	Octadecanol-TMS
2	Adenosine monophosphate-5TMS	44	Glycine-3TMS	86	Octenedioic acid-2TMS
3	Adenosine-4TMS	45	Glycolic acid-2TMS	87	Oleic acid-TMS
4	Alanine-2TMS	46	Glycyl-Glycine-4TMS	88	O-Phosphoethanolamine-4TMS
5	Allose-meto-5TMS(2)	47	Histidine-3TMS	89	O-Phospho-Serine-4TMS
6	4-Aminobutyric acid-2TMS	48	Homocysteine-3TMS	90	Ornithine-3TMS
7	4-Aminobutyric acid-3TMS	49	2-Hydroxyadipic acid-3TMS	91	Ornithine-4TMS
8	3-Aminopropanoic acid-3TMS	50	3-Hydroxybutyric acid-2TMS	92	Oxalic acid-2TMS
9	Arachidonic acid-TMS	51	2-Hydroxyglutaric acid-3TMS	93	5-Oxoproline-2TMS
10	Arginine-3TMS	52	3-Hydroxyglutaric acid-3TMS	94	Palmitic acid-TMS
11	Ascorbic acid-4TMS	53	2-Hydroxyisobutyric acid-2TMS	95	Palmitoleic acid-TMS
12	Asparagine-3TMS	54	3-Hydroxyisobutyric acid-2TMS	96	Pantothenic acid-3TMS
13	Aspartic acid-3TMS	55	3-Hydroxypropionic acid-2TMS	97	Phenylalanine-2TMS
14	Benzoic acid-TMS	56	Hypotaurine-3TMS	98	Phosphoenolpyruvic acid-3TMS
15	Cholesterol-TMS	57	Inositol phosphate-7TMS	99	3-Phosphoglyceric acid-4TMS
16	Citric acid-4TMS	58	Inositol-6TMS(2)	100	Phosphoric acid-3TMS
17	Cystamine-nTMS	59	Isocitric acid-4TMS	101	Proline-2TMS
18	Cystathionine-4TMS	60	2-Isopropylmalic acid-3TMS	102	Putrescine-4TMS
19	Cysteine-3TMS	61	2-Ketoglutaric acid-meto-2TMS	103	Pyridoxine-3TMS
20	Dihydroxyacetone phosphate-meto-3TMS(1)	62	Kynurenine-3TMS	104	Pyruvic acid-meto-TMS
21	Dihydroxyacetone phosphate-meto-3TMS(2)	63	Lactic acid-2TMS	105	Ribitol-5TMS
22	Elaidic acid-TMS	64	Lauric acid-TMS	106	Ribonic acid-5TMS
23	Erythrose-meto-3TMS(1)	65	Leucine-2TMS	107	Ribose-4TMS(2)
24	Fructose 1-phosphate-meto-6TMS(1)	66	Lysine-4TMS	108	Ribulose 5-phosphate-meto-5TMS(1)
25	Fructose 6-phosphate-meto-6TMS	67	Maleic acid-2TMS	109	Saccharopine-4TMS
26	Fructose-meto-5TMS(2)	68	Malic acid-3TMS	110	Serine-3TMS
27	Fumaric acid-2TMS	69	Mannose 6-phosphate-meto-6TMS(1)	111	Sorbitol-6TMS
28	Galacturonic acid-5TMS(2)	70	Mannose 6-phosphate-meto-6TMS(2)	112	Sorbose-meto-5TMS(1)
29	Galacturonic acid-meto-5TMS(1)	71	Mannose-5TMS(1)	113	Stearic acid-TMS
30	Galacturonic acid-meto-5TMS(2)	72	Mannose-meto-5TMS(2)	114	Succinic acid-2TMS
31	Glucaric acid-6TMS	73	Margaric acid-TMS	115	Threonic acid-4TMS
32	Glucose 6-phosphate-meto-6TMS(1)	74	Methionine-2TMS	116	Threonine-3TMS
33	Glucose-5TMS(2)	75	3-Methylglutaconic acid(Z)-2TMS	117	Thymine-2TMS
34	Glucose-meto-5TMS(1)	76	7-Methylguanine-2TMS	118	Trehalose-8TMS
35	Glucose-meto-5TMS(2)	77	Monostearin-2TMS	119	Tyrosine-3TMS
36	Glucuronic acid lactone-3TMS(2)	78	Myristic acid-TMS	120	Urea-2TMS
37	Glucuronic acid-meto-5TMS(1)	79	N-Acetylaspartic acid-2TMS	121	Ureidosuccinic acid-3TMS
38	Glutamic acid-3TMS	80	N-Acetylaspartic acid-3TMS	122	Valine-2TMS
39	Glutamine-3TMS	81	N-Acetylglutamine-3TMS	123	Xylitol-5TMS
40	Glutamine-4TMS	82	N-Acetylmannosamine-meto-4TMS(2)	124	Xylulose-meto-4TMS
41	Glycerol 3-phosphate-4TMS	83	Niacinamide-TMS		
42	Glycerol-3TMS	84	Nonanoic acid-TMS		

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