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
Detection of drugs at different data acquisition speeds by iMScope QT

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



(A) Imipramine, chloroquine, and their possible fragmentation patterns


(B) MS² spectra of standard imipramine

(C) MS² spectra of imipramine detected from the mouse kidneys

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(E) MS² spectra of chloroquine detected from the mouse kidneys

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(A) Imipramine and its metabolites detected in the mouse kidneys

(B) Chloroquine and its metabolites detected in the mouse kidneys

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(A) Imipramine and its metabolites detected in the mouse brains

(B) Chloroquine and its metabolites detected in the mouse brains

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
(A) Rapid detection of imipramine and its metabolites in the mouse kidneys by iMScope QT

(B) Rapid detection of chloroquine and its metabolites in the mouse kidneys by iMScope QT

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<参考文献>

- 1) Buchberger, A.R.; DeLaney, K.; Johnson, J.; Li, L. Mass spectrometry imaging: a review of emerging advancements and future insights. *Analytical chemistry* 2018, 90, 240.
- 2) Harada, T.; Yuba-Kubo, A.; Sugiura, Y.; Zaima, N.; Hayasaka, T.; Goto-Inoue, N.; Wakui, M.; Suematsu, M.; Takeshita, K.; Ogawa, K. Visualization of volatile substances in different organelles with an atmospheric-pressure mass microscope. *Analytical chemistry* 2009, 81, 9153-9157.
- 3) Jackson, S.N.; Muller, L.; Roux, A.; Oktem, B.; Moskovets, E.; Doroshenko, V.M.; Woods, A.S. AP-MALDI mass spectrometry imaging of gangliosides using 2, 6-dihydroxyacetophenone. *Journal of The American Society for Mass Spectrometry* 2018, 29, 1463-1472.
- 4) Doyno, C.; Sobieraj, D.M.; Baker, W.L. Toxicity of chloroquine and hydroxychloroquine following therapeutic use or overdose. *Clinical Toxicology* 2021, 59, 12-23.
- 5) Obuchowicz, E.; Bielecka-Wajdman, A.; Zieliński, M.; Machnik, G.; Gołyszny, M.; Ludyga, T. Imipramine and venlafaxine differentially affect primary glial cultures of prenatally stressed rats. *Frontiers in Pharmacology* 2020, 10, 1687.
- 6) Islam, A.; Takeyama, E.; Mamun, M.A.; Sato, T.; Horikawa, M.; Takahashi, Y.; Kikushima, K.; Setou, M. Green nut oil or DHA supplementation restored decreased distribution levels of DHA containing phosphatidylcholines in the brain of a mouse model of dementia. *Metabolites* 2020, 10, 153.
- 7) Islam, A.; Sakamoto, T.; Zhai, Q.; Rahman, M.M.; Mamun, M.A.; Takahashi, Y.; Kahyo, T.; Setou, M. Application of AP-MALDI Imaging Mass Microscope for the Rapid Mapping of Imipramine, Chloroquine, and Their Metabolites in the Kidney and Brain of Wild-Type Mice. *Pharmaceuticals* 2022, 15, 1314.
- 8) Pazhayattil, G.S.; Shirali, A.C. Drug-induced impairment of renal function. *International journal of nephrology and renovascular disease* 2014, 7, 457.
- 9) Chang, G.-R.; Hou, P.-H.; Wang, C.-M.; Lin, J.-W.; Lin, W.-L.; Lin, T.-C.; Liao, H.-J.; Chan, C.-H.; Wang, Y.-C. Imipramine Accelerates Nonalcoholic Fatty Liver Disease, Renal Impairment, Diabetic Retinopathy, Insulin Resistance, and Urinary Chromium Loss in Obese Mice. *Veterinary sciences* 2021, 8, 189.
- 10) Murugavel, P.; Pari, L. Attenuation of Chloroquine - Induced Renal Damage by α - Lipoic Acid: Possible Antioxidant Mechanism. *Renal failure* 2004, 26, 517-524.
- 11) Yang, J.; Hellerstein, D.J.; Chen, Y.; McGrath, P.J.; Stewart, J.W.; Peterson, B.S.; Wang, Z. Serotonin-norepinephrine reuptake inhibitor antidepressant effects on regional connectivity of the thalamus in persistent depressive disorder: evidence from two randomized, double-blind, placebo-controlled clinical trials. *Brain communications* 2022, 4, fcac100.
- 12) Wichmann, T.O.; Damkier, H.H.; Pedersen, M. A Brief Overview of the Cerebrospinal Fluid System and Its Implications for Brain and Spinal Cord Diseases. *Frontiers in Human Neuroscience* 2021, 15.
- 13) Matsumae, M.; Sato, O.; Hirayama, A.; Hayashi, N.; Takizawa, K.; Atsumi, H.; Sorimachi, T. Research into the physiology of cerebrospinal fluid reaches a new horizon: intimate exchange between cerebrospinal fluid and interstitial fluid may contribute to maintenance of homeostasis in the central nervous system. *Neurologia medico-chirurgica* 2016, 56, 416-441.
- 14) Senova, S.; Fomenko, A.; Gondard, E.; Lozano, A.M. Anatomy and function of the fornix in the context of its potential as a therapeutic target. *Journal of Neurology, Neurosurgery & Psychiatry* 2020, 91, 547-559.

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