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
- 1) Ishimoto S, Fukusaki E, Shimma S. Mass spectrometry imaging of gamma-aminobutyric acid and glutamic acid decarboxylase reactions at various stages of banana ripening, *Journal of Bioscience and Bioengineering* (in press)
- 2) Parijadi AAR, Yamamoto K, Ikram MMM, et al. Metabolome Analysis of Banana (*Musa acuminata*) Treated With Chitosan Coating and Low Temperature Reveals Different Mechanisms Modulating Delayed Ripening, *Frontiers in Sustainable Food Systems*, 6 (2022).
- 3) Marriott J and Palmer JK. BANANAS - PHYSIOLOGY AND BIOCHEMISTRY OF STORAGE AND RIPENING FOR OPTIMUM QUALITY, *Crc Critical Reviews in Food Science and Nutrition*, 13, 41-88 (1980)
- 4) Maduwanthi SDT and Marapana R. Induced Ripening Agents and Their Effect on Fruit Quality of Banana, *International Journal of Food Science*, 2019 (2019)
- 5) Li W, Shao YZ, Chen WX, et al. The Effects of Harvest Maturity on Storage Quality and Sucrose-Metabolizing Enzymes During Banana Ripening, *Food and Bioprocess Technology*, 4, 1273-1280 (2011)
- 6) Wang HB, Zhang ZQ, Xu LY, et al. The effect of delay between heat treatment and cold storage on alleviation of chilling injury in banana fruit, *Journal of the Science of Food and Agriculture*, 92, 2624-2629 (2012)
- 7) Wang YS, Luo ZS, Huang XD, et al. Effect of exogenous  $\gamma$ -aminobutyric acid (GABA) treatment on chilling injury and antioxidant capacity in banana peel, *Scientia Horticulturae*, 168, 132-137 (2014)
- 8) Li Z, Peng Y, and Huang BR. Alteration of Transcripts of Stress-Protective Genes and Transcriptional Factors by  $\gamma$ -Aminobutyric Acid (GABA) Associated with Improved Heat and Drought Tolerance in Creeping Bentgrass (*Agrostis stolonifera*), *International Journal of Molecular Sciences*, 19 (2018)
- 9) Ma Z, Richard H, Tucker DL, Conway T, et al. Collaborative regulation of *Escherichia coli* glutamate-dependent acid resistance by two AraC-like regulators, GadX and GadW (YhiW), *Journal of Bacteriology*, 184, 7001-7012 (2002)
- 10) Salvatierra A, Pimentel P, Almada R, et al. Exogenous GABA application transiently improves the tolerance to root hypoxia on a sensitive genotype of *Prunus* rootstock, *Environmental and Experimental Botany*, 125, 52-66 (2016)
- 11) Xu B, Sai N, and Gilliam M. The emerging role of GABA as a transport regulator and physiological signal, *Plant Physiology*, 187, 2005-2016 (2021)
- 12) Podlesáková K, Ugena L, Spíchal L, et al. Phytohormones and polyamines regulate plant stress responses by altering GABA pathway, *New Biotechnology*, 48, 53-65 (2019)
- 13) Fait A, Nesi AN, Angelovici R, et al. Targeted Enhancement of Glutamate-to- $\gamma$ -Aminobutyrate Conversion in Arabidopsis Seeds Affects Carbon-Nitrogen Balance and Storage Reserves in a Development-Dependent Manner, *Plant Physiology*, 157, 1026-1042 (2011)
- 14) Jin J, Hua L, Zhu ZF, et al. *GAD1* Encodes a Secreted Peptide That Regulates Grain Number, Grain Length, and Awn Development in Rice Domestication, *Plant Cell*, 28, 2453-2463 (2016)
- 15) Inoue K, Shirai T, Ochiai H, et al. Blood-pressure-lowering effect of a novel fermented milk containing gamma-aminobutyric acid (GABA) in mild hypertensives, *European Journal of Clinical Nutrition*, 57, 490-495 (2003)
- 16) Gramazio P, Takayama M, and Ezura H. Challenges and Prospects of New Plant Breeding Techniques for GABA Improvement in Crops: Tomato as an Example, *Frontiers in Plant Science*, 11 (2020)
- 17) Nishimura M, Yoshida S-I, Haramoto M, et al. Effects of white rice containing enriched gamma-aminobutyric acid on blood pressure, *Journal of traditional and complementary medicine*, 6, 66-71 (2016). (in Japanese)
- 18) Akihiro T, Koike S, Tani R, et al. Biochemical mechanism on GABA accumulation during fruit development in tomato, *Plant and Cell Physiology*, 49, 1378-1389 (2008)
- 19) Gallego PP, Whotton L, Picton S, et al. A ROLE FOR GLUTAMATE-DECARBOXYLASE DURING TOMATO RIPENING - THE CHARACTERIZATION OF A CDNA-ENCODING A PUTATIVE GLUTAMATE-DECARBOXYLASE WITH A CALMODULIN-BINDING SITE, *Plant Molecular Biology*, 27, 1143-1151 (1995).
- 20) Liu X, Hu XM, Jin LF, et al. Identification and transcript analysis of two glutamate decarboxylase genes, *CsGAD1* and *CsGAD2*, reveal the strong relationship between *CsGAD1* and citrate utilization in citrus fruit, *Molecular Biology Reports*, 41, 6253-6262 (2014)
- 21) Hu W, Liu JH, Yang XY, et al. Identification of a Gene Encoding Glutamate Decarboxylase Involved in the Postharvest Fruit Ripening Process in Banana, *Hortscience*, 49, 1056-1060 (2014)
- 22) Shelp BJ, Bozzo GG, Trobacher CP, et al. Hypothesis/review: Contribution of putrescine to 4-aminobutyrate (GABA) production in response to abiotic stress, *Plant Science*, 193, 130-135 (2012)
- 23) Lorente-Mento JM, Valero D, Martínez-Romero D, et al. Preharvest Multiple Applications of GABA Improve Quality Traits and Antioxidant Compounds of Pomegranate Fruit during Storage, *Horticulturae*, 9 (2023)
- 24) Enomoto H. Unique distribution of ellagitannins in ripe strawberry fruit revealed by mass spectrometry imaging, *Current Research in Food Science*, 4, 821-828 (2021)
- 25) Takeo E, Fukusaki E, and Shimma S. Mass Spectrometric Enzyme Histochemistry Method Developed for Visualizing *In Situ* Cholinesterase Activity in *Mus musculus* and *Drosophila melanogaster*, *Analytical Chemistry*, 92, 12379-12386 (2020)
- 26) Takeo E, Sugiura Y, Ohnishi Y, et al. Mass Spectrometric Enzyme Histochemistry for Choline Acetyltransferase Reveals De Novo Acetylcholine Synthesis in Rodent Brain and Spinal Cord, *ACS Chemical Neuroscience*, 12, 2079-2087 (2021)
- 27) Wisman AP, Minami M, Tamada Y, et al. Visualization of dipeptidyl peptidase B enzymatic reaction in *rice koji* using mass spectrometry imaging, *Journal of Bioscience and Bioengineering*, 134, 133-137 (2022)
- 28) Ikuta S, Shinohara N, Fukusaki E, et al. Mass spectrometry imaging enables visualization of the localization of glutamate decarboxylase activity in germinating legume seeds, *Journal of Bioscience and Bioengineering*, 134, 356-361 (2022)
- 29) Wang B, Miao Y, Chen J, et al. Activity and distribution of glutamate decarboxylase in banana peel, *Journal of Chinese Institute of Food Science and Technology*, 14, 211-217 (2014) (in Chinese)
- 30) McLellan MR, Lind LR, and Kime RW. HUE ANGLE DETERMINATIONS AND STATISTICAL ANALYSIS FOR MULTIQUADRANT HUNTER L,a,b DATA, *Journal of Food Quality*, 18, 235- 240 (1995)
- 31) Cordenunsi-Lysenko BR, Nascimento JRO, Castro-Alves VC, et al. The Starch Is (Not) Just Another Brick in the Wall: The Primary Metabolism of Sugars During Banana Ripening, *Front Plant Sci*, 10, 391 (2019).

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