FTIR TALK LETTER

vol.9

August 2007

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Jianhua Song¹⁾ and Yuki Hashi²⁾

Shimadzu International Trading (Shanghai) Co., Limited (SSL),

1) Analytical Instrument Division, Marketing Department 2) Customer Support Center

FTIR is one of the major analytical instruments in China nowadays. Of course, the main purpose to use an FTIR instrument is qualitative analysis of organic compounds. It is assumed that the market size of FTIR in China is US\$25M as an annual sales amount and 850 units as a number in year 2006. Last 5 years, FTIR market in China had been expanded, and we are expecting to continue the expansion of this market for another few years. With the growth of FTIR market, apparently, the developments of much more complicated applications are needed as the market requirements. The market segmentation for Shimadzu FTIR in 2006 is shown in Fig.1.

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
Fig.1 FTIR segments in China in 2006

Pharmaceutical industry is the biggest market, followed by University. Because of large population in China, pharmaceutical industry is one of most important industries. Thus, a lot of FTIR instruments have been sold for quality controls of the medicines. Main applications for pharmaceutical industry are identification of drugs as well as its package materials. ATR accessories, which can avoid complex sample pretreatment procedures, are often used to analyze package materials. The package materials are measured with an ATR accessory and their spectra are shown in Fig.2.

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
Fig.2 FTIR spectra of package materials

As the results of search, these materials are identified as polypropylene(Red line in Fig.2) and polyurethane(Blue line in Fig.2). In addition, as a historical reason, there is a specific application for an FTIR instrument in China market. This is the analysis of Traditional Chinese Medicine (TCM). Main ingredients of TCM such as ginsengs may depend on a harvest year, a growth area, etc. Therefore, identification of TCM is getting more complex. However, the use of infrared spectrum makes this identification much easier. Currently, 1103 numbers of infrared spectra are registered in Chinese Pharmacopeia. Compared to British Pharmacopeia and Japanese Pharmacopeia, the number of infrared spectra registered in Chinese Pharmacopeia is much larger than that of other Pharmacopeia. All raw materials and final medicines for TCM must be identified by FTIR. This is the main reason why more than 50 % of market segment in China is occupied by pharmaceutical industry.

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China has been called a factory of the world. In fact, a lot of electronics companies set up the factories in Southern China. The products manufactured by those factories are exported to Europe, USA, Japan, etc. In order to improve the product reliability and comply with some of regulations such as RoHS directive, screening of harmful compounds became big issues. Screening analysis always needs rapid quantitative information. An FTIR instrument can provide some of judgment information for existence of hazardous compounds in the products. In this point, big impact for FTIR application has come from RoHS (Restriction on Hazardous Substances) directive in the last few years. Brominated flame retardants such as Polybromobiphenyls (PBBs) and Polybromodiphenyl ethers (PBDEs) are banned by RoHS directive. To measure these residual PBBs and PBDEs, FTIR is one of effective method to identify these prohibited compounds in the plastic materials. In this regard, electronics industry is also one of rapid growth markets for FTIR instruments.

In all industries, failure analysis is one of the measurement purposes for an FTIR instrument. As an example of failure analysis, the result of contaminant analysis with an ATR accessory is shown in Fig.3.

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
Fig.3 FTIR Spectrum of contaminant

As the result of search, a spectrum of acrylic compound is similar to the contaminant. The contaminant is guessed to be acrylic glue because the contaminant is adhesive. In the case of very small contaminants, microscope-FTIR is getting more demand for failure analysis.

As a factory of the world, FTIR market in China is steadily expanding in the pharmaceutical industry and some of material science industries. In order to meet versatile requirements, applications using FTIR accessories will take an important role.

SSL promotion activities

1. Seminars

Instrument promotion seminars and application seminars were held with some hot topics, such as RoHS and drug analysis in recent years.

2. User Meeting

To improve our after service, FTIR User Meeting was hold in every two years. User1s papers were collected and famous FTIR experts were invited to give presentations. Through the meeting, the relationship between users and Shimadzu was highly strengthened.

3. Academic Conferences

The brand image of Shimadzu FTIR has been improved by attending academic conferences, such as National Molecular Spectrum Conference, where high quality papers by collaboration with famous universities were presented.

4. Internal Sales Training

In order to improve the sales of FTIR, internal training is held regularly. By training, sales become more confident in the competition.

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