nCoV

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P/N:241-09560-91

Instruction Manual

Please read this document carefully and use it correctly. Please store this manual carefully so that it can be used at any time.

Features

- Detection of the N1 and N2 regions on N gene (Nucleocapsid) of the novel coronavirus (SARS-CoV-2) by 1step Reverse Transcription PCR (RT-PCR).
 - This kit uses primers and probe sequences of "2019 - nCoV _ N1" and "2019 - nCoV _ N2" described in "2019 Novel Coronavirus (2019 - nCoV) Real-time rRT-PCR Panel Primers and Probes" by CDC of U.S.
- From sample treatment to RT-PCR in a single tube without RNA extraction.
- Simply mix sample and sample treatment reagent and heat before RT-PCR.
- Contains internal control DNA (IC) to prevent false negative.

Kit contents

No.	Item	Quantity	Volume
①	Sample Treatment Reagent	1	500 μ L
②	nCoV Reagent A	1	650 μ L
③	nCoV Reagent B	1	650 μ L
④	nCoV Reagent C	1	200 μ L

Usage of one kit : 100 tests

Expiry : On the packaging bag label

Storage temperature : -20°C

Other necessities (equipment and materials)

- 1) Real-time PCR system: FAM, ROX and Cy5 detection
- 2) Heat block (adjustable up to 90°C)
- 3) Micropipette and tip with filter
- 4) Small centrifuge for spin down
- 5) Vortex mixer
- 6) Crushed ice or refrigerant
- 7) Aluminum block for tube cooling
- 8) Reaction tube for real-time PCR
- 9) Microcentrifuge tube for reaction mixture (1.5 mL~2 mL)

Procedure

Precautions for use

- Store this kit in a **freezer (-20°C)**.
- All operations should be performed in an ice-cold condition (on aluminum blocks cooled with coolant such as crushed ice).
- Thaw ①Sample Treatment Reagent, ②nCoV Reagent A, ③ nCoV Reagent B to room temperature, mix well with a vortex mixer for 5 seconds, spin down, and store in ice until use.
- ④**Keep** nCoV Reagent C in the freezer **until use (-20°C)**.

1 Sample treatment

- (1) Stir the virus transport medium^{※1} with a vortex mixer for 5 seconds.
 - ※1 Take appropriate measures to prevent coronavirus infection before handling.
- (2) Add 5 μ L of ① Sample Treatment Reagent and 5 μ L of the sample to the PCR reaction tube, mix with a vortex mixer for 5 seconds, and spin down.
- (3) Perform heat treatment for 5 minutes at 90°C.
- (4) Spin down and put on ice.

2 Preparation of reaction mixture

- (1) Prepare solution of RT-PCR in a preparation tube^{※2}.

[RT-PCR solution for one test]

Multiply the required number of reactions to calculate the usage.

② nCoV Reagent A	6.5 μ L
③ nCoV Reagent B	6.5 μ L
④ nCoV Reagent C	2 μ L
Total	15 μL

※2 : Since nCoV Reagent C is a very small amount, be careful when handling it.

After mixing each reagent, mix well with a vortex mixer for 5 seconds.

- (2) Add 15 μ L of the RT-PCR solution "(1)" into a tube containing 10 μ L of the treated sample.
- (3) Mix well with a vortex mixer for 5 seconds, spin down, set the PCR reaction tube in the real-time PCR system, and start the reaction immediately.

3 RT-PCR conditions

- (1) The conditions for RT-PCR are as follows^{※3}.

[Condition]

- Using Thermo Fisher Scientific QuantStudio 5

Temperature	Time
42°C	10min
	↓
95°C	1min
	↓
[95°C	5sec \Rightarrow 60°C 30sec ^{※4} (Detection)]
	×45 cycle

- Using BioRad CFX96 Touch Deep Well

Temperature	Time
42°C	10min
	↓
95°C	1min
	↓
[95°C	5sec \Rightarrow 60°C 15sec ^{※4} (Detection)]
	×45 cycle

※3: Optimization of RT-PCR amplification conditions may be required depending on the Real-time PCR system used. Please validate the results before use. Use the Real-time PCR system in accordance with its operating instructions.

※4: 60°C steps are metered with Cy5, ROX and FAM fluorescent filters. This step should be 30 seconds including the detection time. Set it to 30 seconds depending on the system you use. For example, BioRad CFX96 Touch Deep Well detects after 15 seconds of set time, and total setting time should be 30sec.

- (2) Check the amplification curve (Detection of Cq value) by Cy5, ROX and FAM fluorescent filters. The real-time PCR system may calculate an incorrect Cq value depending on the shape of the amplification curve. Be sure to check the shape of the amplification curve.

【Judgment】

- Confirm that the increase of the amplification curve of the internal control (Cy5) is within 40 cycles. However, if N1 (ROX) or N2 (FAM) increases quickly, internal control may not increase due to lack of PCR substrate.
- If the ROX amplification curve increases within 45 cycles, N1 is judged positive.
- If the FAM amplification curve increases within 45 cycles, N2 is judged positive.

Judgment	ROX: N1 (≤45 cycles)	FAM: N2 (≤45 cycles)	Cy5: IC (≤40 cycles)
Positive	+	+	+ or -
	+	-	+ or -
	-	+	+ or -
Less than detection sensitivity	-	-	+
RT-PCR Fail	-	-	-

+: Increase in amplification curve

-: No increase in amplification curve

Notes

1. Notes on Reagents

- **This kit is for research use only.** It cannot be used for the purpose of diagnosis or for the procedure.
- Follow the instructions in this manual and in Safety Data Sheets.
- This reagent kit does not guarantee complete detection of coronavirus. It may not be detected depending on the amount of virus, amount of contaminants in the sample, or variant virus strains.

2. Notes on Disposal

- To prevent contamination by amplification products, discard the reaction tube after PCR without opening the lid. Do not autoclave before disposal. DNA is not degraded in an autoclave. Aerosols may generate and cause pollution.
- Please dispose of the waste appropriately in accordance with the laws and regulations of the local government.

3. Others

- This manual is copyrighted by Shimadzu Corporation. No part or all of the contents may be reproduced or copied without the permission of Shimadzu Corp.
- The company and product names described in this manual may be trademarks of their respective companies. In this text "®" and "™" is omitted.

About Warranty

1. Warranty

- The expiration date is written on the label of the packaging bag. Please use before the expiration date.
- If there is any problem with the quality before the expiration

date, we will replace the product free of charge.

- The expiration date is valid when the package is unopened. Use promptly after opening.

2. limitation of liability

- In no event shall our company be liable for your lost profits, indirect, or consequential damages. Our company will not be liable for any damages caused to you by third parties.
- In any event, our company's liability for damages shall be limited to the amount equivalent to the cost of the Reagent Product.

3. Exclusions from Warranty

- Abnormalities in quality within the expiration date are excluded from the warranty in the following cases.

(1) If misused

(2) Measurement results become inaccurate due to the customer's operation method.

(3) If the storage method was not appropriate

(4) Abnormalities caused by reasons not related to this product.

Importer in EU:

Shimadzu Europa GmbH
Albert-Hahn-Strasse 6-1047269 Duisburg, F.R. Germany
TEL +49(203)7687-0

Importer in North America:

Shimadzu Scientific Instruments, Inc.
7102 Riverwood Drive Columbia, MD 21046 U.S.A.
TEL +1-800-477-1227

Importer in South and South East Asia:

Shimadzu (Asia Pacific) Pte Ltd.
79 Science Park Drive #02-01/08, Cintech IV
Science Park I Singapore - 118264
TEL (+65) 6778 6280

Manufacturer:

Shimadzu Corporation
1, Nishinokyo Kuwabara-cho, Nakagyo-ku, Kyoto 604-8511, Japan
<https://www.shimadzu.com/>