High Performance Packed Column for HPLC

CoreFocus

Shim-pack PREP Series

INSTRUCTION MANUAL

1. Introductions

- Shim-pack "PREP" series of columns are 20 mm I.D. x 250 mm long stainless steel tubes packed with totally porous, irregular-shaped silica particles (particle size: 15 μ m, pore size: 10 nm) with chemically modified surfaces. This series has been designed exclusively for preparative LC.
- The silica packing is fully end-capped to suppress residual silanol-group influences.

2. Operating Care

- Mobile phase flow in the direction indicated on the column label only.
- Maximum pressure is 25 MPa. Lower pressures (up to 15 MPa) result in a longer lifetime at high efficiency.
- To protect the chemically modified silica surface, do not use solutions of lower than pH 2 or higher than pH 7.5.
- Take care not to apply shocks to the column, such as dropping it on the floor.
- It is recommended to filter the mobile phase and sample solution with a membrane filter or similar beforehand. Contaminants in the mobile phase or sample solution may lead to clogging of the filter at the column's top or clogging of the packing material, resulting in a reduced lifespan of the column.

3. Selection of Mobile Phase

<PREP-ODS, -C8, -TMS>

- A polar organic solvent (e.g., methanol, acetonitrile, etc.) and an aqueous buffer are commonly used as a mobile phase.
- For the separation of ionic substances, appropriate adjustment of pH and/or addition of salts and/or counter ions may result in improved elution times and peak shapes. Recommended are as follows:

Salt: Potassium biphosphate

Counter ion (acidic samples): tetrabutylammonium ions

Counter ion (basic samples): pentanesulphonic ions

Note: Make sure mobile phase pH is within the specified usable range.

- Non-polar substances that do not elute with only methanol or acetonitrile may be eluted by adding tetrahydrofuran or chloroform to methanol or acetonitrile.

<PREP-NH₂>

Reversed Phase Applications

Please use it in the same way as PREP-ODS.

Normal Phase Applications

N-hexane, chloroform, etc. are used, and isopropanol, ethanol, etc. are added to speed up the elution of the compounds.

At the time of shipment, an aqueous solvent is enclosed, so when using it in normal phase, please wash it with isopropanol once. (At this time, be sure to set the flow rate to 2 mL/min or less.)

Anion Exchange Applications

Use a mixture of phosphate buffer, acetonitrile, and methanol.

<PREP-SIL>

N-hexane, chloroform, etc. are used, and isopropanol, ethanol, etc. are added to stop the elution of components.

Column	P/N	Functional group	Guaranteed minimum theoretical plates	Shipping solvent
Shim-pack PREP-SIL	228-00814-91	-	4000	n-Hexane/ethanol = 90/10
Shim-pack PREP-ODS	228-00815-91	Octadecyl group	5000	85% Methanol
Shim-pack PREP-C8	228-00816-91	Octyl group	5000	70% Methanol
Shim-pack PREP-TMS	228-00817-91	Trimethyl group	5000	60% Methanol
Shim-pack PREP-NH ₂	228-17879-91	Aminopropyl group	4000	65% Acetonitrile

4. Maintenance of Column

<General Precautions>

- When removing the column for storage, always seal both ends and store it in a place with minimal temperature fluctuations.
- If the column pressure increases, first remove the column and check for any resistance up to the column inlet. An increase in the column's own resistance is often caused by a blockage in the inlet filter (both the inlet and outlet filters are 2 µm).
- When washing the column, it is recommended to set the flow rate to 2 mL/min or lower.

<PREP-ODS, -C8, -TMS>

- After use, wash and store with methanol.
- If salts have been added to the mobile phase, thoroughly rinse the salts with water first, then store using the above method.
- If amines such as triethylamine or tetra-n-butylammonium have been added to the mobile phase, wash with a methanol/0.5% phosphoric acid solution mixture (50/50, v/v) and then with methanol. (Do not rinse with water only.)
- To remove non-polar substances strongly retained by the packing material, wash with highly soluble solvents such as tetrahydrofuran or chloroform.

<PREP-CN, -NH₂>

In the case of reversed phase mode

Same as in the case of PREP-ODS, etc.

In the case of normal phase mode

- Ensure to enclose and store the mobile phase.
- To eliminate polar compounds that are strongly held by the packing, wash using isopropanol or similar solvents.
- For sugar analysis, wash with a 50/50 (v/v) mixture of acetonitrile and water, and subsequently store in a 70/30 (v/v) mixture of acetonitrile and water.

In the case of anion exchange mode

Flush a water/acetonitrile mixed solvent with the same organic solvent ratio as the mobile phase to wash out the salts. Then store it in a 70/30 (v/v) mixture of acetonitrile and water.

<PREP-SIL>


- Store with the same solvent as the mobile phase.
- For removing adsorbed polar compounds from the column, wash with isopropanol, etc.

5. Column Connection

- For installing the column with stainless steel piping, the following connection parts are required.

Product	P/N	Qty.
Malenut 1.6 MN	228-16001	2
Ferrule 1.6 F	228-16000	2

- Connect the column as illustrated in the figure below. Be careful not to tighten it too much.

6. Disposal Precautions

Separate the column from general industrial or household waste and dispose of it in accordance with local government regulations.

7. Technical Support

If you encounter any performance defects, please contact your local Shimadzu sales office for assistance. However, please note that we cannot guarantee the lifetime of the columns, nor can we accept any claims when performance has degraded due to non-compliance with the operation procedures outlined above, or as a result of normal aging.

*Please note that the appearance, specifications, and the contents of this instruction sheet are subject to change without notice.

High Performance Packed Column for HPLC

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