## Shimadzu Packed Column for HPLC

**CoreFocus**

# Shim-pack™ SUR-Na

## INSTRUCTION MANUAL

### Introduction

Shim-pack SUR-Na is a size-exclusion/ligand-exchange chromatography column used to separate oligosaccharides up to three sugars and monosaccharides.

### Specification

#### Packing Material

Item	Contents
Base material	Sulfonated styrene-divinylbenzene resin
Degree of bridge construction	8%
Particle size	8 µm
Counter ion	Sodium (Na <sup>+</sup> )

#### Analytical Column

P/N	Description	Dimension
228-59529-01	Shim-pack SUR-Na	250 mm × 7.8 mm <i>i.d.</i>

#### Guard Column

P/N	Description	Dimension
228-59529-02	Shim-pack SUR-Na (G)	50 mm × 7.8 mm <i>i.d.</i>

#### Operational Condition

Condition	Contents
Operating pressure	Max. 4.0 MPa
Operating temperature	Max. 85 °C (usually 60-80°C)
Flow rate	Max. 0.8 mL/min (usually 0.6 mL/min)

#### Shipping Solvent : 0.1% Sodium nitrate solution

### Column Performance

The column is shipped after a strict quality control inspection. The results are summarized on the inspection results sheet, which is included with the column(Guard column is excluded).

The performance report includes the column serial number, column performance, etc. Please keep it.

At the time of shipment, the product meets the test specification described in the certificate, but this specification does not guarantee retention time or peak shape for all compounds applicable to this column.

**NOTE:** Compound retention times and peak shape vary with usage. Before developing an analysis method with the columns in this series, verify that the columns satisfy the above standards.

### Column Installation

- The flow direction of the column is shown on the column(→). When installing the column, ensure that this flow direction matches the mobile phase flow direction.
- Use either stainless steel (SUS) or PEEK tubing with an inner diameter of 0.25 - 0.3 mm (0.1 - 0.13 mm for the 4.6 mm column) and an outer diameter of 1.6 mm. Prevent peak broadening caused by dead volume by keeping tubing lengths as short as possible. However, the length of the tubing between the injector and the column needs to be close to 300 mm to allow for preheating.
- To connect the column to the flow lines, use the Male nut, PEEK. Ensure that the fittings are connected properly to avoid creating dead volume between the tubing and the column interface. Extra nuts can be ordered by referring to the part number below.

P/N	P/N	Comment
228-18565-84	Male nut, PEEK	5/pkg


Fig. 1 : Column Connection

- As shown below, the guard column is installed between the injector and the column. Its flow direction arrow also matches the mobile phase flow direction. Use PEEK nuts (male) for the connections.


Fig. 2 : Installing the Guard Column

### Mobile Phase Solvent

- The column is used with water as mobile phase.
- Aqueous solutions of sodium salts can be used.
- Acetonitrile can be mixed into the mobile phase at concentrations below 30%.
- Solvents should be HPLC grade or an equivalent, and reagents should be special reagent grade or an equivalent.

**NOTE:** Never use a mobile phase other than water or the aqueous solutions described.

### Column Handling Precautions

- Observe the usage conditions described in the "Specifications". Set the flow rate and temperature so that the pressure does not exceed 4 MPa.
- Sudden pressure changes lead to deterioration of column performance, so keep the flow rate low until the column temperature rises.
- When the mobile phase is water, sucrose-containing samples analyzed at high temperatures using a SUR-Na column may hydrolyze sucrose and degrade the peak shape. In this case, inject 100 mL of 1 mol/L sodium acetate solution 5 ~ 10 times. If it still does not improve, clean the column according to the procedure described in the "Washing the column".
- Filter the mobile phase, rinse solvent, and sample through a membrane filter (0.2 - 0.45 µm) before use. Suspended particles can clog the column.
- Do not shock the column.

### Washing the column

The adhesion of cations or oily compounds on the column adversely affects both the retention time and peak shape. In such cases, washing the column following the procedure below may recover. Wash the guard column and the analysis column separately.

- Set column temperature at 50°C. Connect the outpipe of the column to the waste bottle. (Do not allow the cleaning solution to flow to the detector side.)
- Flow 0.05 mol/L sodium hydroxide aqueous solution at a flow rate of 0.5 mL/min for 2 hours.
- Flow water at a flow rate of 0.5 mL/min for 2 hours. Ensure that the pH of the solution eluted from the column is neutral and that the sodium hydroxide has been sufficiently washed out.

### Column Storage

When not using the column for a while, remove it from the device, plug it at both ends and store to prevent the inside of the column from drying out.

When storing for longer than a week, store in a cool and dark place and fill the column with the mobile phase. At least once a month, flow the mobile phase at a flow rate of 0.5 mL/min for approximately one hour.

### Technical Support

Shimadzu offers technical support for customers who need to help. Contact your local representative or a nearby Shimadzu sales office if needed.

※ The contents of this instruction sheet are subject to change without notice.

**Shimadzu Corporation**

Analytical & Measuring Instruments Division  
1, Nishinokyo Kuwabara-cho, Nakagyo-ku, Kyoto, 604-8511, Japan

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