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## Instruction Manual

### ■ Introduction

This product is a Shim-pack Fast-OA which is an organic acid analysis column by ion exclusion chromatography. Please read the following instructions before use.

### ■ Specifications

The product specifications are as follows.

Item	Shim-pack Fast-OA(G)
Base material	styrene divinyl benzene copolymer
Ion exchange group	sulfo group
Particle size (μm)	5
Column size	100 mm L. × 7.8 mm I.D.
Max. Flow rate (mL/min)	0.8
Max. Pressure	25.0 MPa or less and 7.0 MPa or less as the differential pressure per 1 column
Number of connectable columns	Guard column and up to 3 analytical column
Max. Temperature (°C)	50 °C

### ■ Operating Precautions

Please check if anything is missing or damaged before use. If there are any signs of damage, notify your local Shimadzu representative at once. The Shim-pack Fast-OA series columns are shipped with the solvent used for the final QC test of the column, as detailed in the column performance report delivered with the column (Guard column and cartridges are excluded). The performance report includes the column serial number, column performance, etc. Please keep it.

At the time of shipment, the product meets the test specification described in the certificate, but this specification does not guarantee retention time or peak shape for all compounds applicable to this column.

**Note** The elution time and peak shape of the compound change depending on the condition of use. Make sure that the column meets the specifications when develop and validate an analytical method using this column.

### ■ Column Installation

The flow direction of the column is shown on the column (→). When installing the column, ensure that the flow direction matches the mobile phase flow direction.

The presence of voids in the tubing connect part may cause leakage and deterioration of column performance (theoretical plate, peak symmetry). Precaution the ferrule tip length or cut surface of the tubing to avoid voids.

Use the shortest tubing connection from the injector to the column to minimize peak broadening.

The column should be connected with male nuts. Male nuts can be ordered by referring to the part number below.

Item name	P/N	Remarks	Pressure
Male nut, PEEK	228-18565-84	5 pcs	20 MPa
UHPLC Fitting 2.5	228-56867-41	1 pc	130 MPa

**Note** Contaminations or air may shorten the column life. Before connecting the column, be sure to flush the flow line by mobile phase.

### ■ Mobile phase

This column basically uses an acid solution as the mobile phase. For example, phosphoric acid, p-toluenesulfonic acid, perchloric acid are used. Keep the following points when preparing the mobile phase.

- Use a reagent with a special grade, HPLC grade or equivalent grade or higher. Use water with HPLC grade or higher. If the metal ion concentration in the reagent is high, separation performance may be degraded. If impurities in the reagent affect the analysis, it may be possible to improve by changing the reagent.

### ■ Sample

When preparing the sample to be injected into this column, pay attention to the following.

- The sample that can be injected into this column is an aqueous solution with pH 1-11. It may contain an water soluble organic solvent, as long as the concentration is 10% or less. If it exceeds these ranges, dilute with water or mobile phase solution.
- If the sample contains highly hydrophobic compounds or compounds that precipitate or gel in the mobile phase, remove the compounds by pretreatment such as solid phase extraction or liquid phase extraction.
- If the sample contains an ionic polymer such as protein, remove it similarly by extraction or ultrafiltration.
- Filter the sample with a membrane filter (0.2 to 0.45 μm).

**Note** Contamination from the sample may accumulate, and may show a increase in column pressure, changes in retention time of the target compound, or deterioration of peak shape. Be sure to connect a guard column when using it and replace it periodically.

### ■ Flushing the Column

Flush the column with 0.1% aqueous phosphoric acid solution at a flow rate of 0.6 mL/min and a column temperature of 50 °C for at least 3 hours. Flush the column for 10 minutes with water, and replace it with the mobile phase.

### ■ Column Handling Precautions

Do not drop or hit the column. Strong shocks can cause the column to deterioration.

To maximize column life, use the columns within the usage condition shown in the “■ Specifications”.

Also, use the P.Max around continuously or suddenly change the pressure, which may cause a decrease in the column life, should be a precaution. Column should be disconnected from the system after the pressure gauge indicates “0”.

### ■ Column storage

When the column will not be used for a while, remove the column from the LC system, plug both ends with a stop plug, and store in a room with minimize temperature fluctuation. Do not remove the cartridge. Do not store in places where temperatures may rise or where there is a risk of freezing. There is no problem that column filling solution keeps the mobile phase when stored for less than one month. When storing for more than one month, replace the filling solution periodically with the mobile phase in the column performance report.

### ■ Disposal Precautions

When disposing of the column, do so in accordance with the processing standards determined by law, separately from general industrial waste and household garbage.

### ■ Technical Support

Should you find any defect in performance, please contact your local Shimadzu representative, who will ensure your complete satisfaction. We regret that we cannot guarantee the lifetime of columns, also that we cannot accept any claim when performance has deteriorated due to noncompliance with the operation procedures elucidated above, or as a result of normal aging.