Shimadzu Packed Column for HPLC

Shim-pack™ SCR-101H

Instruction Manual

■ PRINCIPLE FOR SEPARATION

AND PROCEDURE OF APPLICATION

● Column Packing Material:

This column is packed with a cation exchange resin which is sulfonated poly(styrene-co-divinylbenzene).

● Analyzable Samples:

The COLUMN is available specially for analyses of various organic acids. As some basic materials will not be eluted due to adsorption of them into the COLUMN in some cases, take care not to inject them into the COLUMN unreasonably. It is recommended to use the COLUMN for analyzing various organic acids comprised in fruits, juices, Japanese "SAKE", wines and beverages, etc.

● Mobile Phase:

The COLUMN is filled with an aqueous acidic solution in advance. Use distilled or deionized water adjusted of pH to be acidic as mobile phase. In case the COLUMN is attached to a System used with another mobile phase before, it is necessary to wash the flow line thoroughly and replace it with water completely in advance. After that, wash the flow line with an aqueous acidic solution and then connect the COLUMN to the System.

NOTE: As carelessly letting outside air into the COLUMN to generate bubbles there may cause deterioration of the performance of COLUMN, it is recommended to attach the COLUMN to a System after eliminating any air through the flow line of system.

● Pressure and Flow Rate:

The adequate flow rate of mobile phase for the COLUMN is approximately 0.5 to 1 mL/min. If the flow rate becomes too high, the COLUMN inlet pressure will become too high, resulting in poor separation. If the flow rate becomes too low, the longer separation time will be required. In order to maintain the durability of the COLUMN for a long time, be sure to use it at the flow rate of 0.5 to 1 mL/min mentioned above.

NOTE : The flow rate of 0.5 mL/min will make the pressure of approximately 1 MPa; on the other hand, 1 mL/min will make approximately 2 MPa, varying it more or less depending on the operating temperature of COLUMN and on each COLUMN.

● Temperature:

Although better separation and higher sample loading capacity can be obtained at higher operating temperature with the COLUMN generally, it is recommended to keep operating temperature below 60°C for prevention of the Column from deterioration.

■ MOUNTING AND DISMOUNTING
OF COLUMN ON AND FROM
SYSTEM


Wash the flow line of System well with clean solvent Prior to mounting the Shim-pack SCR-101H on the System. After that, replace the solvent with aqueous acidic solution completely.

As the proper FLOW direction of mobile phase is marked on the COLUMN, never to reverse flow direction through the Column.

Tubing kit available for column connection is listed below.

Item Name	P/N	Comments
Ferrule 1.6F 316L	228-16000-10	1/pkg
Male nut 1.6MN	228-16001	1/pkg

Observe the appropriate connection method as illustrated in the Figure.

Do not overtighten as this will result in damage to the fittings.

Never to change the pressure abruptly. It is required to change the pressure gradually.

Take care not to allow the mobile phase in the reservoir to be used up during the operation for analysis.

Prior to stopping the operation of System, make sure than all the components of the sample injected into the COLUMN have been eluated and the base line of the chromatogram has returned to the initial position. If the analysis has been made at the temperature of COLUMN above room temperature, it is necessary to feed the mobile phase continuously in the same condition as that during the analysis until the temperature of COLUMN has returned to room temperature.

In case of making the analysis successively tomorrow, leave the System as it is.

In case the analysis will not be made for coming 2 to 3 days, disconnect the connection tube between the outlet of the COLUMN and inlet of the detector and then close up the outlet of COLUMN with the Plug 1.6p to prevent the packing material in the COLUMN from drying. In case the COLUMN will not be used for coming fairly long period, it is required to feed the mobile phase through it for about one hour at least once every 4 weeks.


In case the dismounting the COLUMN from the System, close up the both ends of it with the stop plug immediately after its removing to prevent the packing material from drying.

Store the COLUMN in a place which is free from direct sunshine and temperature changes. In cold season, it is recommended to store the column indoors, taking care not to let the packings material freeze.

CAUTION: In case of stopping the operation of System, replace the solution in the reservoir with distilled water and then open The drain valve to clean the flow line of pumping system thoroughly.

NOTE: At least 50 mL of water should be fed to clean the flow line.

In case the analysis will not be made for coming 2 to 3 days, be sure to remove the COLUMN from the System and then clean the flow line throughout the whole system including pump and detector with water thoroughly.

Typical Chromatogram showing the Performance of COLUMN

■ TECHNICAL SUPPORT

Shim-pack SCR-101H columns are manufactured, inspected, packaged and shipped under strict standards of quality control. Should you find any defect in performance, please contact your local Shimadzu representative, who will ensure your complete satisfaction.

We regret that we cannot guarantee the lifetime of columns, also that we cannot accept any claim when performance has deteriorated due to noncompliance with the operation procedures elucidated above, or as a result of normal aging.

※ The contents of this instruction sheet are subject to change without notice.

Shimadzu Corporation

Analytical & Measuring Instruments Division
1, Nishinokyo Kuwabara-cho, Nakagyo-ku, Kyoto, 604-8511, Japan

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