Amino シリーズ

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Test Kit for High-performance Affinity Chromatography

Shim-pack

Amino Series

Instruction Manual

■ Introduction

The Shim-pack Amino Series columns are packed with a strong acid cation exchange resin and suitable for use in the separation of amino acids. Since the column packing is a sulfonic acid ion exchange resin of a styrene-divinyl benzene copolymer (gel type) of 5 μm particle size as a base, this column can be used as a general purpose cation exchange column (if analytical conditions are properly selected).

The Shim-pack Amino Series is classified as follows according to the counter ion of the sulfonic acid group:

Name of Column	Column Size	Counter Ion	Primary Application
Shim-pack Amino-Na	6 mmφ×100 mm	Na type	For analyzing amino acids in protein hydrolysates
Shim-pack Amino-Li	6 mmφ×100 mm	Li type	For separation of amino acids of biological origin

To ensure the best performance and effective use of these columns, please read this operation manual carefully prior to use.

■ Package Content

After unpacking, check that the following accessories are supplied together with the column main body:

- Column Test Record
This record is evidence of the actual performance of the column you have purchased.
- Column connector kit
Containing two male nuts 1.6 MN and two ferrules 1.6 F which are needed for connecting the Shim-pack Amino series column.
- Instruction Manual This manual
(For amino acid analysis, refer to the AMINO ACID ANALYSIS SYSTEM INSTRUCTION MANUAL.)

■ Column Performance

Each Shim-pack Amino Series column is subjected to severe performance tests and only those passing these tests are shipped.

<Amino acid analysis kit>

Kit Name	P/N	Bottle Name	Content	Quantity	Application
Amino acid mobile phase kit NA type	228-21195-94	AA-MA (NA)	0.2N-sodium citrate	1 L	Liquid A, buffer solution No.1
		AA-MB (NA)	0.6N-sodium citrate	1 L	Liquid B, buffer solution No.2
		AA-MC (NA)	0.2N-sodium hydroxide	500 mL	Liquid C, column regenerating solution
Amino acid mobile phase NA type liquid A (3 bottles)	228-21195-96	AA-MA (NA)	0.2N-sodium citrate	1 L × 3 bottles	Liquid A, buffer solution No.1
		AA-MA (LI)	0.15N-lithium citrate	1 L	Liquid A, buffer solution No.1
		AA-MB (LI)	0.3N-lithium citrate	1 L	Liquid B, buffer solution No.2
Amino acid mobile phase kit LI type	228-21195-95	AA-MC (LI)	0.2N-lithium hydroxide	500 mL	Liquid C, column regenerating solution
		AA-MA (LI)	0.15N-lithium citrate	1 L × 3 bottles	Liquid A, buffer solution No.1
		AA-RA	Alkali buffer solution	500 mL	Reagent A
Analyzing kit OPA reagent	228-21195-93	AA-RB	o-Phthalaldehyde solution	500 mL	Reagent B

● Amino-Li type column

Flush the column with 0.2 M lithium hydroxide aqueous solution or mobile phase liquid C (Li type).
(Flush with the cleaning solution at 0.2 mL/min for at least 30 minutes, although these cleaning condition depends on the degree of contamination of the column.)
After the cleaning, run with 0.15 M lithium citrate aqueous solution or mobile phase liquid A (Li type) at 0.6 mL/min for 120 minutes. This cleaning process is necessary to obtain reproducibility of the peak area of sarcosine (Sar).

<Other cleaning method>

● Cleaning with acid

Allow 0.1 N nitric acid to flow through the column at 0.2 mL/min for 1 to 2 hours, substitute the Amino-Na type column with 0.2 M sodium hydroxide and the Amino-Li type column with 0.2 M lithium hydroxide for about few hours and then start analysis.

● Cleaning with EDTA

Allow 0.2 M EDTA-2Na (ethylenediaminetetraacetic acid disodium salt) aqueous solution to flow through the column at 0.2 mL/min for 1 to 2 hours, flush with water, substitute the Amino-Li type column with 0.2 M lithium hydroxide for a day and then start analysis.

If column performance cannot be recovered by the above procedures, it is necessary to replace the column with a new one.

For the column test conditions and actual column performance, please refer to the supplied Test Record.

When shipped, the column is filled with the following solution:
0.2 M sodium hydroxide aqueous solution for the Shim-pack Amino-Na column
0.2 M lithium hydroxide aqueous solution for the Shim-pack Amino-Li column

■ Column Installation

<Precautions for column connection>

Connect the column as shown in Fig.1 using the supplied column connector kit.


Fig.1 : Column Connection

At this time, do not use unnecessarily long pipes for connection between the injector and column and between the column and detector. Excessive pipe length will cause components to disperse widely outside the column (band broadening) so that the column will not appear to provide the expected performance. Standard sizes for pipes are 0.2 to 0.3 mmφ in diameter and 5 to 20 cm in length. (For pipe size for amino acid analysis, refer to the manual for the amino acid analysis system.) Do not tighten male nuts excessively. Excessive tightening could damage joints.

<Direction of mobile phase flow in the column>

Pass mobile phase through the column in the correct direction. Direction of mobile phase flow is marked on the column.

<Flow rate of mobile phase>

The optimum mobile phase flow rate for analysis is between 0.3 and 0.8 mL/min.

<Column operating pressure>

Set the flow rate, ensuring that pressure applied on the column, taking into account the viscosity of the mobile phase, does not exceed 130 kgf/cm².

<Column operating temperature>

Use the column at the temperature not higher than 95 °C.

<Other precautions>

- When analysis is over, drop the mobile phase flow rate to 0.1 to 0.3 mL/min and cool the column to room temperature prior to stopping mobile phase flow.
- Prior to opening the pump drain valve, stop supplying mobile phase and make sure that the column pressure has dropped sufficiently.
- Handle columns carefully so as to avoid shocks (such as from dropping on the floor).

<Analysis conditions for amino acids>

For specific analysis conditions for amino acids, refer to the instruction manual for the amino acid analysis system issued by us.

■ Mobile Phase Solvent

<For amino acid analysis>

In separating amino acids, accurate and proper preparation of mobile phase is essential to obtain reproducible data. To ensure long column life and highly reproducible data, use the appropriate amino acid analysis kit manufactured by us.

<Amino acid analysis kit>

Shimadzu can supply various amino acid analysis kits as shown in the table.

<For separation of substances other than amino acids>

When using the Shim-pack Amino series as a general-purpose cation exchange column to separate substances other than amino acids, take care of the following points.

- Applicable solvent
An ordinary buffer solution can be used without impairing stability. When an organic solvent is added to the mobile phase, select a solvent which is miscible with water. (EX. methanol or acetonitrile) The solvent concentration in the water must not exceed 10 %. Higher concentrations might cause contraction of the packing gel.
- pH range
pH of the mobile phase must be in the range between 2 and 14.
- Filtration of mobile phase
To ensure long, stable performance of the column, pass the mobile phase through a membrane filter (of 0.2 to 0.45 μm pore size) prior to use.

■ Sample

To prevent a clogged column, filtrate the sample solution through a membrane filter (0.2 - 0.45 μm), if possible.

■ Column Handling Precautions

Observe the following to assure long service life of the column:

- Rinse the flow line in the LC system sufficiently prior to connecting the column.
If soil accumulating in the flow line in the LC system flows in the column, it will deteriorate column performance. Therefore, be sure to rinse the flow path and piping system of the pump and injector sufficiently before connecting the column. For rinsing, use an organic solvent such as methanol, isopropyl alcohol or distilled water according to the operating conditions of the LC system.

<For amino acid analysis>

Arrange the system according to the applications guide for the amino acid analysis system issued by us. After connecting the piping, clean the flow line.

■ Clogging of Column

Increase in column pressure during analysis may be caused by any of the following:

- (1) Flow line before the column inlet is clogged.
- (2) Sintered filter in the column inlet is clogged.
- (3) Column packing bed is clogged.
- (4) Flow line after the column outlet is clogged.

If column pressure increases, disconnect the pipe from the column inlet and check for item (1).

If item (1) is not the cause, reconnect the pipe to the column inlet and disconnect the pipe from the column outlet to check for item (4). If a clogged flow line is not the cause, either item (2) or (3) is the cause of the column pressure rise.

When item (3) is the cause, the problem may be eliminated by cleaning the column in the procedure described in "■ Cleaning of the Column".

When item (2) is the cause, disconnect the column end and subject the filter embedded in the column end to ultrasonic cleaning. Then the problem may be eliminated. In disconnecting the column end, take care not to disturb the packing bed; if the bed is disturbed, column performance will deteriorate. (The filter is embedded in the column end. Therefore, perform ultrasonic cleaning of the filter as embedded in the column end.)

■ Flushing the Column

Irreversible adsorption of sample components on the column packing may result in retention time change and/or poor resolution. In such a case, the following column cleaning procedure may be effective.

<When amino acid analysis has been conducted:>

- Amino-Na type column
Flush the column with 0.2 M sodium hydroxide aqueous solution or mobile phase liquid C (Na type).
(Flush with the cleaning solution at 0.2 mL/min for at least 30 minutes, although these cleaning condition depends on the degree of contamination of the column.)
After the cleaning, run with mobile phase liquid A at 0.5 mL/min for 20 minutes.