# Shimadzu High Performance Packed Column for HPLC

**Shim-pack**

## IC-SA3/-SA3(G)

### Instruction Manual

#### ■ Introduction

Shim-pack IC-SA3 is a polymer gel-based column with an anion exchange group for ion chromatography. It is a high-resolution analytical column used for the analyses of anions such as chlorate, chlorite, bromate in addition to fluoride, chloride, bromide, nitrate, nitrite, sulfate and phosphate. This column is applicable to suppressed ion chromatography. Shim-pack IC-SA3(G) is a guard column for protection of Shim-pack IC-SA3.

#### ■ Specifications

##### ● Packing

Item	Contents
Type of support	Polyvinyl alcohol resin
Ion exchanger	Quaternary ammonium group
Particle size	5 $\mu$ m

##### ● Construction Materials

Material
PEEK (polyether ether ketone)

##### ● Analytical Column

P/N	Description	Dimensions
228-41600-91	Shim-pack IC-SA3	250 mm $\times$ 4.0 mm <i>i.d.</i>

##### ● Guard Column

P/N	Description	Dimensions
228-41600-92	Shim-pack IC-SA3(G)	10 mm $\times$ 4.6 mm <i>i.d.</i>

##### ● Operational Conditions

Condition	Contents
Pressure	Max. 15 MPa
Temperature	20 - 60 °C
Flow rate	Max. 1.0 mL/min
pH	3 - 12
Organic solvent	Max. 10% either of methanol or acetonitrile

When adding an organic solvent to the mobile phase, check the specifications of the equipment including the suppressor. Note that organic solvents cannot be used with the cartridge suppressor and membrane suppressor (ICDS-40A) manufactured by Shimadzu (as of July 2023).

#### ■ Column Performance

Each Shim-pack IC-SA3 column is shipped only after strict factory inspection. The results are summarized on the performance report with the test chromatogram, which is included in the package. (It is not applicable to a guard column.) Shim-pack IC-SA3 is supposed to meet the right standard and performance as indicated below, but this does not represent a guarantee of the retention time or peak shape for every applicable samples.

**NOTE:** Compound retention times and peak shape may vary with usage. Before developing an analytical method with the column, verify that the column satisfies the above standards.

case of the sample containing silicate, a negative peak associated with silicate may be observed before fluoride. In addition, it is possible to adjust retention of carbonate and/or silicate by optimizing the mobile phase condition.

Description	Criteria	Chromatographic Conditions
IC-SA3	$N > 10,000$ ( $\text{NO}_2^-$ )	Mobile Phase: 3.6 mmol/L Sodium Carbonate Flow Rate: 0.8 mL/min Temperature: 45 °C Detection: Suppressed Electroconductivity Sample: Nitrite ion 5 $\mu$ g/mL Injection Volume: 50 $\mu$ L

#### ■ Column Installation

Consider the following points when this series of the column is installed on the system.

- The flow direction of the column is shown on the column (→). When installing the column, ensure that this flow direction matches the mobile phase flow direction.
- Use PEEK tubing with an inner diameter of 0.25 - 0.3 mm and an outer diameter of 1.6 mm. Prevent peak broadening caused by dead volume by keeping tubing lengths as short as possible. However, the pre-heating coil kit (P/N 228-45714- 91) must be set between the injector and the column.
- The column is connected with supplied PEEK Male nuts. Ensure that the fittings are connected properly to avoid creating dead volume between the tubing and the column interface. Extra nuts can be ordered by referring to the part number below.
- The guard column is installed between the injector and the analytical column. Use supplied PEEK Male nuts for the connections.

Item Name	P/N	Comment
Male nut, PEEK	228-18565-84	5/pkg

**NOTE:** The stain or air in the flow line may deteriorate the column. Before connecting the column, be sure to flow the mobile phase to flush the flow line.

#### ■ Mobile Phase Solvent

This column is generally used with sodium carbonate buffer as a mobile phase. Typical mobile phase is as follows.

3.6 mmol/L Sodium Carbonate aqueous solution.

The higher the pH and/or the concentration of carbonate ions in the mobile phase, the earlier the anions elute. Therefore, it is possible to adjust retention of anions by changing pH or carbonate concentration, if necessary. The pH range of the mobile phase should be 3-12. The mobile phase may contain up to 10 % methanol or acetonitrile. In addition, solvents should be HPLC grade or an equivalent, and the reagents should be high purity grade or an equivalent.

## ■ Sample

The sample which can be injected is an aqueous solution with pH 3-12. Water miscible organic solvents can be included in the sample, only if the concentration is less than 10 %. When the sample does not meet the requirement, it should be diluted with water or mobile phase, or performed a sample pretreatment by a solid-phase extraction. In case of samples containing organic solvents, pressure may increase. For this case, it should be diluted with water or mobile phase as well to prevent from the excessive pressure beyond the maximum limitation.

When the sample is solid which cannot be solved in water, or which is an organic solvent which is non-miscible with water, sample pretreatment must be performed to extract target ions to water. Consider the following points when injecting the sample.

- When the sample contains highly hydrophobic compounds or compounds that precipitate or gel in the mobile phase, remove them by a pretreatment such as solid phase extraction or liquid phase extraction.
- When ionic macromolecules such as proteins are included, remove them from the sample prior to analysis by extraction or ultrafiltration.
- Filter the sample through a membrane filter (0.2 - 0.45 µm) prior to injection.

**NOTE:** The filters specified to ion chromatography should be used. General use filters may contain anions and it may cause contaminate the samples.

## ■ Column Handling Precautions

Consider the following points when using the column.

- Observe the pressure, temperature, and flow rate limits given in “■ Specifications”. In particular, do not allow the pressure to exceed 70 %. The steep pressure change over the column may cause deterioration.
- Disconnect the column at room temperature without applying pressure.
- Do not shock the column by banging it or dropping it.

## ■ Flushing the Column

Some kinds of problems may be occurred such as increase of the column pressure, variation of the retention time, deterioration of peak shape and reduction of the peak area, when metal ions (including hydroxide formation etc.), multivalent anions such as proteins or hydrophobic organic compounds are adsorbed in the column. In case that these phenomena are observed, it might be ameliorable by flushing the column.

- Flushing the column is performed by the appropriate solvent selected one of the three solvents described below based on the primary factor of the contamination. Deliver the appropriate flushing solvent onto the column with the pump.
- The method of cleaning the guard column and the analytical column is the same but wash them separately.
- Place the tube from the exit of the column into a waste bottle directly in order to prevent the waste solution from flowing into the detector.
- Set the half of standard analytical flow rate or 0.5 mL/min as the flow rate for column flushing, unless otherwise specified. Column pressure may be higher depending on the contamination or flushing solvent. In case of such higher pressure, decrease the flow rate to reduce pressure under the typical pressure.
- Typical flushing period should be one to two hours. After

flushing, deliver enough of the mobile phase to wash out the flushing solvent completely.

### Flushing solvent -1

Retrailing of F<sup>-</sup> peak or diminution of PO<sub>4</sub><sup>3-</sup> peak intensity.

Deliver mixture of aqueous solution containing 50 mmol/L EDTA-2Na with acetonitrile (10/1, v/v) for approximately 12 hours. (Subsequent to that, deliver the mobile phase for 15 minutes at 0.5 mL/min followed by at the analytical flow rate for around 2 hours. Then, connect to the suppressor.)

### Flushing solvent -2

For hydrophobic (multivalent electrolytic) contaminants

Deliver 10-fold concentration of the mobile phase. (Subsequent to that, deliver the mobile phase for an hour at 0.5 mL/min and then, connect to the suppressor.)

### Flushing solvent -3

For hydrophobic contaminants

Deliver 5% acetonitrile aqueous solution for 10 minutes followed by acetonitrile for an hour. (Observe the column pressure since it may increase. Subsequent to this procedure, deliver purified water for 30 minutes at 0.5 mL/min and then, the mobile phase for an hour at the same flow rate followed by connecting to the suppressor.)

100% acetonitrile is only for column flushing. For normal use, the maximum of acetonitrile is 10%. Column flushing should be performed only when it is necessary. Arbitrary flushing with different type of solvents may cause reduction of column performance or rise in pressure due to degradation of the resin. Further, these column flushing procedures may not guarantee recovery of the column performance.

**NOTE:** For suppressor system, prior to the column flushing process, place the tube from the exit of the column into a waste bottle directly in order to prevent the waste solution from flowing into the suppressor. After the flushing, deliver enough of the mobile phase to wash remaining solution before the suppressor is connected. In case that solution containing EDTA flows into a suppressor, the suppressor may be clogged.

**NOTE:** The column cannot be regenerated if it is heavily contaminated. Use the guard column and replace it at regular interval to prevent from contamination.

## ■ Column Storage

When the column is not used in a short term, disconnect it from the instrument and cap both ends with the provided stop-plugs and keep in a room under stable ambient temperature. Do not store in the place where it may be under higher temperature or the mobile phase may freeze. If column is not used and stored within one month, it is not necessary to replace the mobile phase. In case of long-time storage more than one month, it is recommended to replace the filled mobile phase into new one at regular interval.

**NOTE:** The packing material cannot be regenerated if it is completely dried. Seal both ends of the column firmly for storage dried. Seal both ends of the column firmly for storage.

## ■ Technical Support

It is the customer's responsibility to develop and validate analytical conditions for a particular application. However, Shimadzu offers technical support by e-mail and phone for customers who need help.

Write specific questions to [analytic@group.shimadzu.co.jp](mailto:analytic@group.shimadzu.co.jp) or call your local representative.

\* The contents of this instruction sheet are subject to change without notice.

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