High Performance Packed Column  
for High Performance Liquid Chromatograph

Shim-pack

IC-A3/-GA3

## Instruction Manual

## ■ Introduction

Shim-pack IC-A3 is an polymer gel based column with an anion exchange group for ion chromatography. It is optimized for non-suppressed ion chromatography with high theoretical plate number and excellent selectivity, and enables the separation of many kinds of organic and inorganic anions. Shim-pack IC-GA3 is a guard column for protection of Shim-pack IC-A3.

## ■ Specifications

## ● Packing

Item	Contents
Type of support	Polyacrylate resin
Ion exchanger	Quaternary ammonium group
Particle size	5 μm

## ● Construction Materials

Material
PEEK (polyether etherketone)

## ● Analytical Column

P/N	Description	Dimensions
228-31706-91	Shim-pack IC-A3	150 × 4.0 mm <i>i.d.</i>

## ● Guard Column

P/N	Description	Dimensions
228-31706-92	Shim-pack IC-GA3	10 × 4.6 mm <i>i.d.</i>

## ● Operational Conditions

Condition	Contents
Pressure	Max. 12.0 MPa (less than 8.0 MPa is recommended)
Temperature	Max. 50 °C (40 °C is recommended)
Flow rate	Max. 1.8 mL/min (1.2 mL/min is recommended)
pH	3 - 8.5
Organic solvent	Max. 70% either of ethanol, methanol or acetonitrile

## ■ Column Performance

Each Shim-pack IC-A3 column is shipped only after strict factory inspection. The results are summarized on the performance report with the test chromatogram, which is included in the package. (It is not applicable to a guard column.) Shim-pack IC-A3 is supposed to meet the right standard and performance as indicated below, but this does not represent a guarantee of the retention time or peak shape for every applicable samples.

 **NOTE:** Compound retention times and peak shape may vary with usage. Before developing an analytical method with the column, verify that the column satisfies the above standards.

Description	Criteria	Chromatographic Conditions
IC-A3	$N > 5,000$ (NO <sub>3</sub> <sup>-</sup> peak)	Mobile Phase: 8 mmol/L <i>p</i> -Hydroxy benzoic acid 3.2 mmol/L Bis-Tris 50 mmol/L Boric acid Flow Rate: 1.2 mL/min Temperature: 40 °C Detection: Non-suppressed Electroconductivity Sample: Nitrate ion 30 μg/mL Injection Volume: 10 μL injection

## ■ Column Installation

Consider the following points when this series of the column is installed on the system.

- The flow direction of the column is shown on the column (→). When installing the column, ensure that this flow direction matches the mobile phase flow direction.
- Use PEEK tubing with an inner diameter of 0.25 - 0.3 mm and an outer diameter of 1.6 mm. Prevent peak broadening caused by dead volume by keeping tubing lengths as short as possible. However, the pre-heating coil kit (P/N 228-45714-91) must be set between the injector and the column.
- The column is connected with supplied PEEK Male nuts. Ensure that the fittings are connected properly to avoid creating dead volume between the tubing and the column interface. Extra nuts can be ordered by referring to the part number below.
- The guard column is installed between the injector and the analytical column. Use supplied PEEK Male nuts for the connections.

Item Name	P/N	Comment
Male nut, PEEK	228-18565-84	5/pkg

## ■ Mobile Phase Solvent

For an anion analysis by a non-suppressed ion chromatography, aromatic carboxylic acid (phthalic acid, *p*-hydroxy benzoic acid etc.) or boric acid can be used as an mobile phase. Typical weak acidic mobile phase condition is shown as below.

- Mixture solution of 8 mmol/L *p*-hydroxy benzoic acid, 3.2 mmol/L Bis-Tris and 50 mmol/L boric acid.

The pH range of the mobile phase for this column should be 3-8.5. It is possible to adjust retention of anions by optimizing concentration of the carboxylic acid in the mobile phase, if necessary. The mobile phase may contain up to 70% either of ethanol, methanol or acetonitrile. In addition, solvents should be HPLC grade or an equivalent, and the reagents should be high purity grade or an equivalent. When impurities in the reagent influence the analysis, it might be ameliorable by changing the regent.

## ■ Sample

The sample which can be injected onto the column is an aqueous solution with pH 3-8.5. Water miscible organic solvents can be included in the sample, only if the concentration is less than 70 %. When the sample does not meet the requirement, it should be diluted with water or mobile phase, or performed a sample pretreatment by a solid-phase extraction. In case of samples containing organic solvents, pressure may increase. For this case, it should be diluted with water or mobile phase as well to prevent from the excessive pressure beyond the maximum limitation. When the sample is solid which can not be solved in water, or is an organic solvent which is non-miscible with water, sample pretreatment must be performed to extract target ions to water. Consider the following points when injecting the sample.

- When highly hydrophobic compounds or compounds with which precipitation or gelation takes place in the mobile phase is included, sample pretreatment such as solid-phase or liquid-phase extraction is necessary prior to analysis.
- When ionic macromolecules such as proteins are included, remove them from the sample prior to analysis by extraction or ultrafiltration.
- Filter the sample through a membrane filter (0.2 - 0.45 μm) prior to injection.

 **NOTE:** The filters specified to ion chromatography should be used. General use filters may contain anions and it may cause contamination in the sample.

 **NOTE:** First time of using column  
This column filled with 70% ethanol after the inspection is shipped. Pressure may be higher with same flow rate of standard mobile phase, as viscosity of ethanol is higher than standard mobile phase. Please keep the flow rate under 1/4 of typical flow rate until they filled solvent in the column have been substituted to the mobile phase to prevent from the over-pressure in the column. It may take 30 min at the flow rate 0.2 mL/min.

## ■ Column Handling Precautions

Consider the following points when using the column.

- Observe the conditions indicated in “■ Specifications”. In particular, do not allow the pressure to exceed 70 %. The steep pressure change over the column may cause deterioration.
- Disconnect the column at ambient temperature without applying pressure.
- Do not shock the column by banging it or dropping it.

## ■ Flushing the Column

Some kinds of problems may be occurred such as increase of the column pressure, variation of the retention time, deterioration of peak shape and reduction of the peak area, when metal ions (including hydroxide formation etc.), multivalent anions such as proteins or hydrophobic organic compounds are adsorbed in the column. In case that these phenomenon are observed, it might be ameliorable by flushing the column.

- Flushing the column is performed by the appropriate solvent selected one of the three solvents described below based on the primary factor of the contamination. Deliver the appropriate flusing solvent onto the column with the pump.

- Place the tube from the exit of the column into a waste bottle directly in order to prevent the waste solution from flowing into the detector.
- Set the half of standard analytical flow rate or 0.5 mL/min as the flow rate for column flushing, unless otherwise specified. Column pressure may higher depending on the contamination or flushing solvent. In case of such higher pressure, decrease the flow rate to reduce pressure under the typical pressure.
- Typical flushing period should be one to two hours. After flushing, deliver enough of the mobile phase to wash out the flushing solvent completely.

<Flushing solvent -1: For multivalent ionic contaminants>  
120 mol/L potassium hydrogen phthalate in aqueous solution

<Flushing solvent -2: For hydrophobic contaminants>  
Mixture of aqueous solution containing 8 mmol/L *p*-hydroxy benzoic acid and 3.2 mmol/L Bis-Tris, with acetonitrile (1/1, v/v) Methanol can be used as a substitute for acetonitrile.

<Flushing solvent-3: For higher pressure or variation of the retention time etc.>  
Mixture of 20 mmol/L (sodium) phosphate buffer pH 2.5 containing 100 mmol/L sodium perchlorate, with ethanol (1/2,v/v). Deliver this solution in the opposite direction against the analytical flow.

Column flushing should be performed only when it is necessary. Arbitrary flushing with different type of solvents may cause reduction of column performance or rise in pressure due to degradation of the resin. Further, these column flushing procedures may not guarantee recovery of the column performance.

 **NOTE:** The column cannot be regenerated if it is heavily contaminated. Use the guard column and replace it at regular interval to prevent from contamination.

## ■ Column Storage

When the column is not used in a short term, disconnect it from the instrument and cap both ends with the provided stop-plugs and keep in a room under stable ambient temperature. Do not store in the place where it may be under higher temperature or the mobile phase may freeze. If column is not used and stored within one month, it is not necessary to replace the mobile phase. In case of long time storage more than one month, it is recommended to replace the filled mobile phase into new one at regular interval.

 **NOTE:** The packing material cannot be regenerated if it is completely dried. Seal both ends of the column firmly for storage.

## ■ Technical Support

It is the customer’s responsibility to develop and validate analytical conditions for a particular application. However, Shimadzu offers technical support by e-mail and phone for customers who need help.  
Write specific questions to [analytic@group.shimadzu.co.jp](mailto:analytic@group.shimadzu.co.jp) or call your local representative.

\* The contents of this instruction sheet are subject to change without notice.

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