High Performance Packed Column for HPLC

CoreFocus

Shim-pack IC-C1/IC-GC1

INSTRUCTION MANUAL

1. Introduction

Shim-pack IC-C1/GC1 is a polymer gel based column for ion chromatography. Shim-pack IC-C1 is a cation exchange column optimized for non-suppressor ion chromatography. By choosing the eluent, it can analyze alkali metal ions, alkylamines, alkaline earth metal ions, and transition metal ions. Shim-pack IC-GC1 is a guard column used to protect the Shim-pack IC-C1.

2. Specification

(Materials)

Packing material	Base material	Styrene-Divinylbenzene copolymer
	Functional group	Sulfo-group
	Particle size	10 μ m
Column body	Material	Stainless steel

(Size)

Parts Number	Product	Size
228-17737-91	Shim-pack IC-C1	150 mm x 5.0 mm I.D.
228-17738-91	Shim-pack IC-GC1	10 mm x 4.0 mm I.D.

(Condition)

Pressure	< 5.0 MPa (Less than 2.0 MPa is recommended)
Temperature	< 50°C (40°C is recommended)
Flow rate	< 20 mL/min (1.5 mL/min is recommended)
pH	2 - 11
Organic solvent	Up to 10% of either Acetonitrile or Methanol can be added.

3. Column Performance

At the time of shipment, each product conforms to the following specifications (excluding guard columns). Please note that these specifications do not guarantee retention time or peak shape for all compounds that can be applied to this column.

Specification	Condition
N > 2,500 (Na ⁺ peak)	Eluent: 5 mmol/L Nitric acid aqueous solution Flow rate: 1.5 mL/min Temp.: 40°C Detection: Conductivity detector (Non suppressor) Sample: Na ⁺ (10 mg/L, 25 μ L injection)

Note Retention times and peak shape may vary with usage. Before developing an analytical method with this column, verify that the column meets the above standard specification.

4. Column Installation

Consider the following points when this column is installed to the LC system.

- The flow direction of the column is shown on the column (→). When installing the column, ensure that this flow direction matches the mobile phase flow direction.
- Use PEEK tubing with an inner diameter of 0.25 - 0.3 mm and an outer diameter of 1.6 mm. To prevent peak broadening caused by dead volume, keep tubing lengths as short as possible. To avoid temperature fluctuations, set the pre-heating coil kit (P/N 228-45714-91) between the injector and the column.
- Use the supplied PEEK male nuts for connecting this column. Ensure that the fittings are connected properly to avoid dead volume between the tubing and the column. Extra nuts can be ordered by referring to the part number below.

Product	Part Number	Remarks
Male nut PEEK	228-18565-84	5 pcs / 1 pack

- The guard column is installed between the injector and the analytical column. Use the supplied PEEK male nuts for the connections.

5. Eluent

The eluent for this column differs according to the valence of the cations in the target ions to be analyzed. The standard eluents are listed below. Please choose one according to the sample.

- Monovalent cations (alkali metals): 5 mmol/L Nitric acid aqueous solution
- Multivalent cations (alkaline earth metals): 4 mmol/L tartaric acid and 2 mmol/L ethylenediamine* mixed aqueous solution

*If the retention is weak, 1 mmol/L ethylenediamine is recommended.

The higher the concentration of hydrogen ions in the eluent, the faster the elution of cations. Therefore, the concentration of acid is adjusted as necessary to control the retention of cations. The eluate used in this column has an acceptable pH range of 2 to 11. The eluent may contain up to 10% methanol or acetonitrile.

- Use reagents of special reagent grade or higher. Specifically, use nitric acid with a low metal ion content for harmful metal analysis. Also, ensure that the water and organic solvents used are of HPLC grade or higher. If impurities in the reagents affect your analysis, consider changing the reagents to improve the results.

Note In the analysis of monovalent cations, install the pre-column (P/N 228-17744-91 Shim-pack IC-PC1) between the pump and the injector to eliminate unnecessary substances (multivalent cations) from the eluent.

6. Analytical Sample

The sample which can be injected onto the column is an aqueous solution with pH2-11. Water miscible organic solvents can be included in the sample, only if the concentration is less than 10%. When the sample does not meet the requirement, it should be diluted with water or mobile phase or performed a sample pretreatment by a solid-phase extraction.

When the sample is solid which can not be solved in water or is an organic solvent which is non-miscible with water, sample pretreatment must be performed to extract target ions to water. Consider the following points when injecting the sample.

- When highly hydrophobic compounds or compounds with which precipitation or gelation takes place in the mobile phase is included, sample pretreatment such as solid-phase or liquid-phase extraction is necessary prior to analysis.
- When ionic macromolecules such as proteins are included, remove them from the sample prior to analysis by extraction or ultrafiltration.
- Filter the sample through a membrane filter (0.2 - 0.45 µm) prior to injection.

Note For sample filtration, please use a membrane filter specifically designed for ion chromatography. Other membrane filters may contaminate the sample with adhered cations.

7. Column Handling Precautions

Consider the following points when using this column.

- Follow the conditions indicated in "Specifications". In particular, the pressure should not exceed 70%. Sudden pressure changes across the column may cause deterioration.
- Disconnect the column at ambient temperature without applying pressure.
- Do not shock the column by banging it or dropping it.

8. Column Cleaning

When hydrophobic substances such as multivalent cations or organic materials in the sample adhere to the column, it may lead to increased column pressure, changes in elution times of target ions, deterioration of peak shape, and reduced peak area values. If such phenomena are observed, check to see if removing the guard column, if it is connected, restores it. In this case, clean or replace only the guard column. If not, clean the analytical column according to the following procedure. The column may recover through cleaning. After cleaning, thoroughly replace the cleaning solution with the eluent to ensure no cleaning solution remains in the flow path.

<Cleaning phase>

20 mmol/L ethylenediamine and 40 mmol/L tartaric acid mixed aqueous solution

- The flow rate during cleaning should be about half the analytical flow rate or 0.5 mL/min. The column pressure may increase depending on the type of contamination and cleaning. In this case, reduce the flow rate to avoid exceeding the normal operating pressure.
- The guard column can be washed in the same way, but the analytical column and the guard column must be washed separately.
- Ensure that the eluent from the column does not flow into the detector. Connecting the piping to a waste bottle directly.
- The cleaning time should be about 1 to 2 hours. After cleaning, thoroughly replace the cleaning solution with the eluent to ensure no cleaning solution remains in the flow path.

Perform column cleaning only when necessary. Frequent cleaning may alter the properties of the packing material, leading to performance degradation and increased pressure. Please note that the cleaning methods described here do not guarantee the restoration of column performance.

Note Depending on the contamination level of the column, the cleaning effect may not be sufficient. Please make sure to always connect a guard column and replace it regularly.

9. Column Storage

When the column is not used, disconnect it from the instrument and cap both ends with the supplied stop-plugs and storage under stable ambient temperature. Do not store in the place where it may be under higher temperature or the eluent may freeze. If column is not used and stored within one month, it is not necessary to replace the analytical eluent. In case of long-time storage more than one month, it is recommended to replace the filled eluent into new one at every month.

Note The packing material cannot be regenerated if it is dried. Ensure that both ends of the column are securely capped during storage.

10. Disposal Precautions

Separate the column from general industrial or household waste and dispose of it in accordance with local government regulations.

11. Technical Support

If you encounter any performance defects, please contact your local Shimadzu sales office for assistance. However, please note that we cannot guarantee the lifetime of the columns, nor can we accept any claims when performance has degraded due to non-compliance with the operation procedures outlined above, or as a result of normal aging.

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