## High Performance Packed Column for HPLC

**CoreFocus**

# Shim-pack IC-A1/IC-GA1

## INSTRUCTION MANUAL

### 1. Introduction

The Shim-pack IC-A1/GA1 is an ion chromatography column filled with a polymer gel that contains anionic exchange groups. The Shim-pack IC-A1 is an analytical column optimized for non-suppressed ion chromatography, designed to analyze inorganic anions and organic acid anions. This column is characterized by its wide usability of eluent pH, ranging from acidic to alkaline conditions. Shim-pack IC-GA1 is a guard column used to protect the Shim-pack IC-A1.

### 2. Specification

(Materials)

Packing material	Base material	Poly(methacrylate) gel
	Functional group	Quaternary ammonium groups
	Particle size	12.5 $\mu\text{m}$
Column body	Material	Stainless steel

(Size)

Parts Number	Product	Size
228-17733-91	Shim-pack IC-A1	100 mm x 4.6 mm I.D.
228-17734-91	Shim-pack IC-GA1	10 mm x 4.6 mm I.D.

(Condition)

Pressure	< 5.0 MPa (Less than 2.0 MPa is recommended)
Temperature	< 50°C (40°C is recommended)
Flow rate	< 20 mL/min (1.5 mL/min is recommended)
pH	2 - 11
Organic solvent	Up to 10% of either Acetonitrile or Methanol can be added.

### 3. Column Performance

At the time of shipment, each product conforms to the following specifications (excluding guard columns). Please note that these specifications do not guarantee retention time or peak shape for all compounds that can be applied to this column.

Specification	Condition
N > 2,000 ( $\text{NO}_3^-$ peak)	Eluent: 2.5 mmol/L Phthalic acid + 2.4 mmol Tris (hydroxymethyl)aminomethane + 16.2 mmol/L Boric acid solution Flow rate: 1.5 mL/min Temp.: 40°C Detection: Conductivity detector (Non suppressor) Sample: $\text{NO}_3^-$ (10 mg/L, 35 $\mu\text{L}$ injection)

**Note** Retention times and peak shape may vary with usage. Before developing an analytical method with this column, verify that the column meets the above standard specification.

### 4. Column Installation

Consider the following points when this column is installed to the LC system.

- The flow direction of the column is shown on the column (→). When installing the column, ensure that this flow direction matches the mobile phase flow direction.
- Use PEEK tubing with an inner diameter of 0.25 - 0.3 mm and an outer diameter of 1.6 mm. To prevent peak broadening caused by dead volume, keep tubing lengths as short as possible. To avoid temperature fluctuations, the pre-heating coil kit (P/N 228-45714-91) is set between the injector and the column.
- Use the supplied PEEK male nuts for connecting this column. Ensure that the fittings are connected properly to avoid dead volume between the tubing and the column. Extra nuts can be ordered by referring to the part number below.

Product	Part Number	Remarks
Male nut PEEK	228-18565-84	5 pcs / 1 pack

- The guard column is installed between the injector and the analytical column. Use the supplied PEEK male nuts for the connections.

### 5. Eluent

In the analysis of anions using non-suppressed ion chromatography, weak acid aqueous solutions such as aromatic carboxylic acids (phthalic acid, p-hydroxybenzoic acid, etc.) and boric acid are generally used as eluents. The standard eluent are as follows.

- 2.5 mmol/L phthalic acid and 2.4 mmol/L tris(hydroxymethyl)aminomethane
- This column can be used with eluents with a pH range of 2 to 11, allowing for analysis under basic eluent conditions. The pH and organic acid concentration in the eluent can be adjusted as needed to modify the retention of anions. Additionally, the eluent may contain either methanol or acetonitrile at concentrations up to 10%.
- Use reagents of special reagent grade or higher. For water and organic solvents, use HPLC grade or equivalent quality. If impurities in the reagents affect the analysis, changing the reagents may improve the results.

### 6. Analytical Sample

The sample which can be injected onto the column is an aqueous solution with pH 2-11. Water miscible organic solvents can be included in the sample, only if the concentration is less than 10%. When the sample does not meet the requirement, it should be diluted with water or mobile phase or performed a sample pretreatment by a solid-phase extraction. When the sample is solid which can not be solved in water or is an organic solvent which is non-miscible with water, sample pretreatment must be performed to extract target ions to water. Consider the following points when injecting the sample.

- When highly hydrophobic compounds or compounds with which precipitation or gelation takes place in the mobile phase is included, sample pretreatment such as solid-phase or liquid-phase extraction is necessary prior to analysis.

- When ionic macromolecules such as proteins are included, remove them from the sample prior to analysis by extraction or ultrafiltration.
- Filter the sample through a membrane filter (0.2 - 0.45µm) prior to injection.

**Note** For sample filtration, please use a membrane filter specifically designed for ion chromatography. Other membrane filters may contaminate the sample with adhered anions.

## 7. Column Handling Precautions

Consider the following points when using this column.

- Follow the conditions indicated in "Specifications". In particular, the pressure should not exceed 70%. Sudden pressure changes across the column may cause deterioration.
- Disconnect the column at ambient temperature without applying pressure.
- Do not shock the column by banging it or dropping it.

## 8. Column Cleaning

Adsorption of metal ions (including hydroxides), multivalent anions such as proteins, and hydrophobic substances from organic materials onto the column can lead to increased column pressure, changes in elution times for target ions, degraded peak shapes, and reduced peak areas. In such cases, performing a column cleaning may partially restore functionality.

- Cleaning should be conducted by selecting one of the two cleaning phase below based on the cause of contamination and flushing the solution using solvent delivery pump.
- The same cleaning method applies to the guard column, but the main column and guard column should be cleaned separately.
- The cleaning phase should be flushed through the column in the same "forward" direction (→) as during analysis.
- Ensure that the eluent from the column does not flow into the detector. Connecting the piping to a waste bottle directly.
- The flow rate during cleaning should be approximately half of the analytical flow rate or 0.5 mL/min. Depending on the column's contamination or the type of cleaning solution, column pressure may increase; in such cases, reduce the flow rate to avoid exceeding the typical operating pressure.
- The cleaning time should be about 1 to 2 hours. After cleaning, thoroughly replace the cleaning solution with the eluent to ensure no cleaning solution remains in the flow path.

### Cleaning Solution 1)

When contaminated by polyvalent electrolytes:  
120 mmol/L potassium hydrogen phthalate aqueous solution.

### Cleaning Solution 2)

When contaminated by hydrophobic substances:  
A mixture of standard eluent and acetonitrile (9/1, v/v)  
(methanol can also be used instead of acetonitrile).

Perform column cleaning only when necessary. Frequent cleaning may alter the properties of the packing material, leading to performance degradation and increased pressure. Please note that the cleaning methods described here do not guarantee the restoration of column performance.

**Note** Depending on the contamination level of the column, the cleaning effect may not be sufficient. Please make sure to always connect a guard column and replace it regularly.

## 9. Column Storage

When the column is not used, disconnect it from the instrument and cap both ends with the supplied stop-plugs and storage under stable ambient temperature. Do not store in the place where it may be under higher temperature or the eluent may freeze. If column is not used and stored within one month, it is not necessary to replace the analytical eluent. In case of long-time storage more than one month, it is recommended to replace the filled eluent into new one at every month.

**Note** The packing material cannot be regenerated if it is dried. Ensure that both ends of the column are securely capped during storage.

## 10. Disposal Precautions

Separate the column from general industrial or household waste and dispose of it in accordance with local government regulations.

## 11. Technical Support

If you encounter any performance defects, please contact your local Shimadzu sales office for assistance. However, please note that we cannot guarantee the lifetime of the columns, nor can we accept any claims when performance has degraded due to non-compliance with the operation procedures outlined above, or as a result of normal aging.

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