High Performance Packed Column for HPLC

CoreFocus

Shim-pack NovaCore C18-HB [Reversed phase column]

INSTRUCTION MANUAL

■ Introduction

This is the instruction manual for the Shim-pack NovaCore C18-HB. The Shim-pack NovaCore C18-HB is a core-shell column that utilizes an organic silica hybrid . To ensure stable and long-term use of the Shim-pack NovaCore C18-HB, which has excellent column performance, please read this instruction manual carefully and use it correctly.

■ Operating Precautions

- Please check for any abnormalities in the appearance and packaging box of the column.
- Verify that there are no errors in the name of the packing material, particle size, column size, etc., by checking the labels attached to the packaging box and the column.
- The inspection report includes the column lot number, column serial number, and chromatographic test conditions. Please keep the report for future reference.
- The number of theoretical plate stated in the inspection report is for quality checks at the time of shipment and does not indicate the performance (specifications) of the column.

■ Column Performance

The Shim-pack NovaCore C18-HB undergoes rigorous QC testing to ensure quality and stability. The shipping solvent is Acetonitrile / Water = 45 : 55 (v : v). Refer to the table below for the specification of the Shim-pack NovaCore C18-HB columns.

Item	Specification
Functional Group	Octadecylsilyl (ODS) group
Pore Size	10 nm
Surface Area	200 m ² /g
Carbon Content	11%
Endcapping	○
Column hardware	Stainless (SUS)
Usable pH range	1-12
Common temperature	20 - 40°C
Maximum temperature	60°C*

* The upper temperature limit is 60°C, but using it at conditions where the pH of mobile phase exceeds 8 at 60°C will accelerate the degradation rate of the column.

NOTE

Near the pH limit of use, please use a mobile phase containing at least 10% organic solvent. At the pH limit, column lifetime may be shortened depending on conditions such as temperature and mobile phase composition. Additionally, continuing to use the column at the maximum temperature limit may also reduce its lifetime.

■ Column Installation

- The flow direction of the column is indicated by an arrow on the column label (→). When installing the column, ensure that the flow direction arrow matches the mobile phase flow direction.
- Use PEEK or SUS tubes with an inner diameter of 0.25 – 0.3 mm (UHPLC: 0.1 – 0.2 mm) and an outer diameter of 1.6 mm. However, do not use PEEK tubes when working with fluorinated organic solvents such as tetrahydrofuran (THF), chloroform, concentrated sulfuric acid, concentrated nitric acid, dichloroacetic acid, acetone, dichloromethane, dimethyl sulfoxide (DMSO), and hexafluoroisopropanol (HFIP).
- The 1.7 μm particle size column produces higher backpressure than the 5 μm or 2.6 μm particle size columns. Please be aware of the maximum pressure of the HPLC system and connect the appropriate tubes. Generally, UHPLC systems that have a maximum pressure above 60 MPa are appropriate for 1.7 μm columns.
- Use the shortest possible tubes connection from the injector to the column to minimize peak broadening.
- Use male nuts or fittings for connecting the column. Be careful to avoid creating any unnecessary gaps during the connection. For information on column connections, please also refer to the system guide and the operating manual of the LC system.
- Male nuts or fittings can be obtained with the following product names and part numbers (P/N).

Item Name	P/N	Remarks	Pressure Limit
Male Nut, PEEK	228-18565-84	5 pcs	20 MPa
Male Nut Fitting Kit	228-45717-01	2 pcs	35 MPa
UHPLC Fitting 2 S	228-56867-41	1 pc	130 MPa

NOTE

Particulates or air in the system flow line may deteriorate the column. Before connecting the column, be sure to filter all solvents and flush the flow line up to the column with mobile phase.

- If peaks are tailing more on the early eluting compounds than later eluting compounds, it's possible that there is dead volume. Make sure the connection tube is fully inserted into the column joint. As a guideline, tighten the connection when the tube extends about 5 mm from the tip of the ferrule, and when removed, the tubes should extend 2-3 mm from the tip of the ferrule.
- Make sure to use appropriate internal diameter and tube length from the injector and to the detector, especially when using semi-micro size columns, to reduce system dead volume and peak broadening.

■ Column Handling Precautions

- To avoid damage to the column and prevent deterioration in performance, do not drop or bump the column.
- To maximize column lifetime, use the column within the pressure limit shown in the following table.

Particle Size	Column I.D.	Maximum Pressure
1.7 μm	1.0 mm	103 MPa*
	2.1 mm	
2.6 μm	2.1 mm	
	3 mm	
	4.6 mm	
5 μm	2.1 mm	
	3 mm	
	4.6 mm	

*Continuous use near the maximum pressure limit can shorten the column's lifetime. Pressure varies depending on factors such as column length, column temperature, and the type of organic solvent, so please adjust the flow rate accordingly.

- The column should be disconnected from the system only after the pressure shows 0 MPa.
- Sudden pressure fluctuations can lead to early degradation of the column. Please avoid slow operation of the sample injection valve.
- If there are issues such as poor reproducibility of compound retention times, baseline drift, or noise at the start of the analysis, it may indicate that the column has not been adequately equilibrated. It is recommended to flush the column with at least five times its volume using the mobile phase intended for the analysis until improvements are observed. If no improvements are seen, please consider reviewing the equipment and analytical conditions.

NOTE

The elution order, retention time, peak shape, etc. may change significantly when ion-pairing reagents are used in the mobile phase. These reagents are often difficult to completely remove, so we recommend that columns with a history of using ion-pair reagents be used exclusively for analyses with ion-pair reagents.

■ Mobile Phase Selection

- Both aqueous and non-aqueous solvents can be used as mobile phases; however, repeatedly switching between solvents with extremely different polarities may lead to a decline in column performance. Common organic solvents that can be used include acetonitrile, methanol, and tetrahydrofuran (THF). However, when using THF, please pay attention to the solvent resistance of the materials used for the tubes and other components in the flow path.

NOTE

When switching from aqueous solvents to non-aqueous solvents, it is recommended to flush with an intermediate polarity solvent, such as isopropanol or dioxane, at least ten times the column volume (approximately 25 mL for a 150 × 4.6 mm I.D. column) before making the substitution. Additionally, since intermediate polarity solvents may result in higher pressure, please be cautious and reduce the flow rate by half during the flushing process to manage the pressure.

- Conditions such as pH and temperature extremes affect the column lifetime. Typical operating temperature is between 20°C and 40°C. When using a column at a high pH for a long time, using low concentration organic buffer solution from 1 to 10 mmol/L at low temperatures is recommended (e.g., < 30°C).

■ Column Flow Rate

The recommended flow rates are as follows.

Particle Size	Column I.D.	Recommended Flow Rate
1.7 μm	1 mm	0.1 mL/min
	2.1 mm	0.5 mL/min
2.6 μm	2.1 mm	
	3 mm	0.8 mL/min
	4.6 mm	1.8 mL/min
5 μm	2.1 mm	0.3 mL/min
	3 mm	0.5 mL/min
	4.6 mm	1.3 mL/min

■ Sample Considerations

Samples should be dissolved in a solvent with a composition as similar as possible to that of the mobile phase (initial solvent during gradient elution). Injecting a large volume of samples dissolved in a solvent with higher solubility than the mobile phase can lead to decreased separation efficiency or precipitation of the sample at the column inlet. Additionally, to prevent the precipitation of salts contained in the sample or the sample solvent within the column, please confirm their miscibility with the mobile phase before injection.

■ Column Clogging

The causes of pressure rise and peak splitting may include clogging or contamination of the filter at the column inlet.

- Filter the mobile phase using a 0.45 μm membrane filter before flushing the flow line and connecting the column.
- Filter samples with a 0.2 – 0.45 μm membrane before injecting.

- Possible causes of baseline drift and noise include pump malfunction due to dissolved air, decreased light intensity when using a UV detector, bubble formation during high-temperature analysis, and solvent purity. If baseline drift and noise occur and do not improve after flushing the column, please consider reviewing the equipment and analytical conditions.
- If macromolecules such as proteins or polysaccharides are adsorbed onto the column, they are generally difficult to remove by cleaning. For samples containing these substances or samples with many impurities, it is recommended to perform pre-treatment in advance.
- During long-term storage, replace the shipping solvent with Acetonitrile/Water = 45:55 (v:v). When storing, seal with the provided plug and keep it in a cool, clean, dark place with minimal temperature fluctuations and low humidity.

■ Precautions When Using Small ID and Small Particle Size UHPLC Columns

The extra-column volume has a significant effect on sample diffusion, especially with 1.0/2.1 mm ID columns. When using this size column, optimize the LC system as described below.

- 1) The tubes from injector - column and column - detector should be as short as possible to minimize dead volume. Recommended tube ID is 0.15 mm or less. Take care when installing the tube to the inlet and outlet ports that no voids are formed in the connection.
- 2) Use a semi-micro or micro flow cell in the UV or PDA detector. Use a fixed volume sample loop to minimize system volume from the injector.
- 3) The detector response and the data sampling rate of the data processing device should be optimized according to the peak width.

■ Cleaning and Storing the Column

In general, columns should be cleaned using the following methods. The flushing volume during cleaning and replacement should be at least five times the column volume.

- If the mobile phase does not contain buffers or salts, increase the concentration of the organic solvent that makes up the mobile phase to wash away substances with significant retention remaining in the column. Organic solvents can be used up to 100%. If highly lipophilic components are adsorbed, adding THF may be effective.
- If the mobile phase contains buffers, salts, or ion-pairing agents, first replace it with a water/organic solvent mixture (in the same ratio as the mobile phase) that does not contain these components, and then clean as described above. For buffers or salts around 50 mmol/L, direct replacement with a 60% acetonitrile aqueous solution is possible. During this process, take care to ensure that salts do not precipitate.
- After using the column near the pH limits, washing with water alone may cause column degradation. It is recommended to clean with the aforementioned water/organic solvent mixture or a 60% acetonitrile aqueous solution.

■ Disposal Precautions

When disposing of the column, do so in accordance with the processing standards determined by law, separately from general industrial waste and household garbage.

■ Technical Support

Shim-pack NovaCore C18-HB columns are manufactured, inspected, packaged and shipped under strict standards of quality control. Should you find any defect in performance, please contact your local Shimadzu representative, who will ensure your complete satisfaction.

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