High Performance Packed Column for HPLC

CoreFocus

Shim-pack™ NovaCore Guard Column [for analysis]

INSTRUCTION MANUAL

■ Introduction

Shim-pack NovaCore Guard Column (for analysis) is a cartridge-type analytical guard column filled with packing material under high pressure. Please read this instruction manual carefully and use it correctly.

■ Operating Precautions

- Please check for any abnormalities in the appearance of the guard column and its packaging.
- Verify that there are no errors in the name of the packing material, particle size, column size, etc., by checking the label attached to the exterior.
- Ensure that there are no scratches or debris on the tip of the cartridge column.

■ Column Performance

The product specifications for the analytical cartridge column are as follows.

Items	P/N	Particle size (μm)	Functional group	Pressure limit (MPa)
Shim-pack NovaCore Guard Column Cartridge, 2.1 mm I.D. x 2 mm, 3 pcs	227-32904-01	5.0	Octadecyl silyl (ODS)	138
Shim-pack NovaCore Guard Column Cartridge, 3.0 mm I.D. x 2 mm, 3 pcs	227-32904-02			
Shim-pack NovaCore Guard Column Cartridge, 4.6 mm I.D. x 2 mm, 3 pcs	227-32904-03			

A separate cartridge holder (P/N: 227-32904-04) is required for installing the cartridge column.

■ Column Installation

1. Confirm that the cartridge holder is composed of two components: the tube base and the floating nut (Fig. 1).

Fig. 1 Component structure of the cartridge holder.

2. Attach the tube base to the floating nut. At this time, ensure that the end of the tube protrudes more than 4 mm from the ferrule (Fig. 2).

Fig.2 Assembly of the tube base and floating nut.

3. Insert the cartridge column into the cartridge holder (Fig. 3) and hand-tighten it clockwise.

Fig.3 Connection direction of the cartridge column.

4. Insert the tip of the tube of the cartridge holder, which is connected to the cartridge column, straight into the end fitting of the analytical column (Fig. 4) and hand-tighten it clockwise.


Fig.4 Connect the cartridge holder to the analytical column.

5. Position the guard column and the analytical column vertically. Insert one wrench into the flat surface of the end fitting of the analytical column and insert another wrench into the flat surface of the cartridge column (Fig. 5). Rotate the wrench on the cartridge column clockwise by approximately 30° to 90°. At this time, do not rotate the wrench inserted into the end fitting of the column. Over-tightening may cause damage.


Fig.5 Tighten the cartridge column with a wrench.

■ Replacement of the Cartridge Column

1. Insert one wrench into the flat surface of the cartridge column and another wrench into the flat surface of the cartridge holder (Fig. 6). Rotate counterclockwise to remove it.


Fig. 6 Removal of the cartridge column.

2. Insert the cartridge column into the cartridge holder and hand-tighten it clockwise.
3. Insert one wrench into the flat surface of the cartridge holder and another wrench into the flat surface of the cartridge column.
4. Rotate the wrench clockwise by approximately 30° to 90°. Over-tightening may cause damage.

■ Usage Precautions

- Repeatedly removing the cartridge column can compromise its airtightness, leading to the risk of liquid leakage.
- The timing for replacing the cartridge column should be determined by monitoring increases in pressure, changes in peak shape, and alterations in separation performance.
- Injecting large amounts of samples repeatedly can cause impurities to adsorb onto the main column, leading to its degradation. It is recommended to replace the cartridge column promptly.

■ Disposal Precautions

When disposing of the column, do so in accordance with the processing standards determined by law, separately from general industrial waste and household garbage.

■ Technical Support

Shim-pack NovaCore Guard Columns are manufactured, inspected, packaged and shipped under strict standards of quality control. Should you find any defects in performance, please contact your local Shimadzu representative, who will ensure your complete satisfaction.

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