Shim-pack

IC-Bromate/(G)

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Shim-pack**IC-Bromate/(G)****Instruction Sheet****1. Introduction**

Shim-pack IC-Bromate/-Bromate(G) is the anion exchange column, which packing material is polymer gel. These are used for separation of bromate ion, chlorous acid ion and iodate ion.

Shim-pack IC-Bromate(G) is the special guard column for Shim-pack IC-Bromate.

2. Specifications

Packing material	Particle	Polyvinyl alcohol gel
	Modification	Quaternary ammonium
	Particle size	9 μ m
Column body	material	PEEK(poly ether ether ketone)

P/N	Product Name	Size
228-46884-91	Shim-pack IC-Bromate	150 mm X4.0 mm i.d.
228-46884-93	Shim-pack IC-Bromate(G)	10 mm X4.6 mm i.d.

Pressure	IC-Bromate	Max7.0 MPa (typ.<3.0 MPa)
	IC-Bromate(G)	Max0.8 MPa (typ.<0.5 MPa)
Temperature		Max60°C (typ.40°C)
Flow rate		Max2.0mL/min (typ.1.0mL/min)

3. Column Performance

The Shim-pack IC-Bromate column is produced under good quality control. The performance report with the test chromatogram is included in this column package (except for the guard column). Shim-pack IC-Bromate/(G) being shipped with the following specification.

Product name	Specification	Condition
IC-Bromate	N> 3500(Bromate) Elution time:5.30~5.70min Column pressure:<3.0MPa	Mobile phase:1.8mM Sodium Carbonate /1.7mM Sodium hydrogen carbonate Flow rate:1.0 mL/min Temperature:40°C Detect:UV(210nm) Sample:Bromate ion, Chlorite ion each 10 μ g/mL Injection volume 50 μ L
IC-Bromate(G)	Column pressure : <0.5MPa	

Note:The elution time and the peak shape vary in the condition. When the development of the analysis method and verification are given by using this column, confirm the column specification before use.

4. Column Installation

When use this column for the bromate analysis system,

please refer to the instruction manual of kit for a bromate analysis system. General information for handling are shown in the following.

- The flow direction of the column is shown on its body (→). When installing the column, ensure that this flow direction matches the mobile phase flow direction.
- Tubings of the connection must be used the PEEK tubing, 0.25 or 0.3 mm I.D. and 1.6 mm O.D.. Do not make them too long to prevent the broadening of the peak. But, only the length of the tubings from the injector to the column should be 300~600mm for preheating.
- The column is connected with supplied PEEK malenuts. Ensure that the fittings are connected properly to avoid creating dead volume between the tubing and the column interface. The product name and the parts number of the spare PEEK malenuts are as follows.

Item name	P/N	Comment
Male nut, PEEK	228-18565-84	5/pkg

- The guard column is installed in front of the column for the analysis.

5. Mobile Phase

This column usually uses sodium carbonate buffer. A typical mobile phase is as following.

1.8mM sodium carbonate / 1.7mM sodium hydrogen carbonate solution.

Anion is eluted earlier if the carbonate ion concentration of the mobile phase is higher, so it is possible to change the elution time of the anion, by pH or the concentration of the sodium carbonate solution adjusting. pH of the mobile phase which can be used for this column is 3 to 12.

Pay attention to the following point when you make a mobile phase.

- Water and an organic solvent to be used must be a HPLC grade. Or, use the thing of the grade according to that.

6. Samples

Samples which can inject into this column are the water solution and pH range is 3 to 12. Dilute it with water or mobile phase, or else process it with the ion removal solid phase extraction column when the sample exceed these ranges.

Pay attention to the following point when you inject a sample.

- If the sample includes highly hydrophobic substances or substances which will be precipitation or gel in the

mobile phase, remove them using by solid phase extraction or dilute by water or the mobile phase.

• Samples are preferably filtered through a membrane filter of 0.2-0.45µm pore size.

7. Handling

Pay attention to the following point when use this column.

- This column should be used under the condition mentioned in the specifications. Especially, please use under 70% of the conditions of the maximum pressure. Please pay attention not to make a rapid pressure change. It will make a column deterioration.
- When remove this column, confirm that there is no pressure and the temperature of the column is as same as room temperature.
- Do not make a shock to the column dropping, ex.

NOTE: The column cannot be regenerated, if the column packings dry. It must be sealed tightly with the supplied plugs when storage.

8. Flushing the column

If the column adsorbs metal ion (including hydroxide), multivalent anion or hydrophobic, the phenomenon such as rise of the column pressure, unstable retention time, bad shape of the peak and reduction of the area value will appear. Flush the guard column if the phenomenon is recovered by removing the guard column from the flow line when such a phenomenon appears.

If the phenomenon is not recovered by removing the guard column from the flow line, please flush both of the analytical column and the guard column. To flush the column, reverse the direction of flow and flush the column as following.

When flush the column, disconnect the column from the flow line, and drain flushing liquid (the drain tubing can be used).

After flushing the column, reconnects the column so as to be the same flow direction as the mobile phase flow direction. To prevent flushing liquid flows in the flow line, replace it with a mobile phase.

If the column do not recover by flushing procedure 1, try flushing procedure 2.

< Flushing procedure 1>

Flow rate	: 0.5mL/min(reverse direction of flow)
5min	: Pure water
60min	: 10mM EDTA-2Na aqueous
5min	: Pure water
Flow rate	: 1.0mL/min(normal direction of flow)
120min	: Mobile phase

< Flushing procedure 2>

Flow rate	: 0.5mL/min(reverse direction of flow)
5min	: Pure water
60min	: Mobile phase of ten times concentration/ acetonitrile (HPLC grade) = 7:3 (v/v)
5min	: Pure water
Flow rate	: 1.0mL/min(normal direction of flow)
60min	: Mobile phase

NOTE: The column cannot be regenerated, if it is heavily contaminated. Use the guard column to prevent from contamination.

9. Column Storage

For storage more than one week, remove the column from the system, cap both ends of the column so that the solvent can not evaporate and keep in the place where a temperature change is small. Do not keep it in the place where temperature in high. When storage the column, mobile phase can be used for sealing liquid. A sealing liquid in shipping is mentioned in [3.Column Preformance].

NOTE: The column cannot be regenerated, if the column packings dry. It must be sealed tightly with the supplied plugs when storage.

10. Technical Support

It is the customer's responsibility to develop and validate analytical conditions for a particular application. However, Shimadzu offers technical support by e-mail for customers who need help.

Write specific questions to

analytic@group.shimadzu.co.jp

in English or Japanese, or call your local representative.

The contents of this instruction sheet are subject to change without notice.

SHIMADZU CORPORATION

**ANALYTICAL & MEASURING
INSTRUMENTS DIVISION**

1, Nishinokyo-Kuwabaracho, Nakagyo-ku, Kyoto, 604-8511, Japan