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The screenshot displays the IMAGEREVEAL software interface. A central dialog box titled "解析手法の選択" (Select Analysis Method) is open, listing several analysis options. The "画像解析" (Image Analysis) option is highlighted with a yellow border and a mouse cursor. The background interface includes a sidebar with various tool icons, a main workspace with a dark image, and a right-hand panel with input fields for parameters like "m/z±許容幅" and "化合物名/コメント".

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The screenshot displays the IMAGEREVEAL software interface. The left sidebar contains a vertical menu of icons. The top icon, labeled "IMDXファイルを追加" (Add .imdx file), is highlighted with a yellow box and a hand cursor. Below it are icons for "画像設定" (Image Settings), "画像レジストレーション" (Image Registration), "ROI設定" (ROI Settings), "まとめて解析" (Batch Analysis), "データ行列" (Data Matrix), "前処理設定" (Preprocessing Settings), "前処理" (Preprocessing), "データ行列設定" (Data Matrix Settings), "データ行列計算" (Data Matrix Calculation), "画像解析" (Image Analysis), "画像分類計算" (Image Classification Calculation), "画像分類結果表示" (Image Classification Result Display), and "類似画像抽出計算" (Similar Image Extraction Calculation).

The main workspace is divided into several panels:

- ROIリスト** (ROI List): A table with columns "No.", "使用" (Use), "ファイル名" (File Name), "ROI名" (ROI Name), and "データ点数" (Data Points).
- データ行列テーブル** (Data Matrix Table): A table with columns "No.", "使用" (Use), "タグ" (Tag), "ラベル" (Label), "m/z", "組成式" (Chemical Formula), "アダクトイオン" (Adduct Ion), "マトリックス" (Matrix), and "極性" (Polarity).
- MSイメージ** (MS Image): A panel showing a dark image with a control bar at the top and a metadata panel on the right with fields for "m/z±許容幅" (m/z±Tolerance), "化合物名/コメント" (Compound Name/Comment), "アダクトイオン" (Adduct Ion), "ファイル名" (File Name), and "タイプ" (Type).
- グラフ** (Graph): A panel with a control bar and a table with columns "計.." (Calculation), "表.." (Table), and "ファイル名" (File Name).
- MSイメージ一覧** (MS Image List): A panel with a control bar and a table with columns "計.." (Calculation), "表.." (Table), and "ファイル名" (File Name).

The "解析パラメータ" (Analysis Parameters) panel at the bottom left shows a table with columns "No.", "Name", and "Value", and a sub-section for "ノーマライズ未計算" (Not Normalized).

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
The screenshot displays the software's main interface. On the left, the 'File' menu is open, and the 'Reference Image Management' option is highlighted with a yellow box and a mouse cursor. The main workspace is divided into several panels: a large empty area at the top, a mass spectrum graph in the bottom-left, and a heatmap image in the top-right. The mass spectrum graph shows intensity versus m/z with several peaks labeled. The heatmap image shows a color-coded distribution of data points. On the right side, there is a panel with metadata for the selected file, including the file name and type.

No.	Name	Value

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

The screenshot displays the JAGEREVEAL software interface. On the left sidebar, the '画像レジストレーション' (Image Registration) option is highlighted with a yellow box and a mouse cursor. The main window is divided into several panels:

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The screenshot displays the IMAGEREVEAL software interface. The sidebar on the left contains several icons, with the '前処理設定' (Pre-processing Settings) icon highlighted by a yellow box and a mouse cursor. The main workspace is divided into several panels:

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The '前処理設定' (Pre-processing Settings) window is currently open, showing a '解析パラメータ' (Analysis Parameters) section with a 'ノーマライズ未計算' (Normalization not calculated) button.

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

The screenshot displays the IMAGEREVEAL software interface with the following components:

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The screenshot shows the IMAGEREVEAL software interface with the following components:

- Left Sidebar:** A vertical menu with various analysis options. The option "類似画像抽出計算" (Similar Image Extraction Calculation) is highlighted with a yellow box and a hand cursor.
- ROIリスト (ROI List):** A table listing ROI information.
- データ行列テーブル (Data Matrix Table):** A table showing data points for each ROI.
- MSイメージ (MS Image):** A large color-coded mass image of the sample.
- グラフ (Graph):** A mass spectrum plot showing intensity versus m/z.
- MSイメージ一覧 (MS Image List):** A list of MS images with a thumbnail.
- Bottom Panel:** A table of analysis parameters.

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The screenshot displays the IMAGEREVEAL software interface with several panels:

- ROIリスト (ROI List):** A table listing ROI parameters.
- データ行列テーブル (Data Matrix Table):** A table listing detected peaks with their m/z values and PLS coefficients.
- MSイメージ (MS Image):** A large color-coded image showing the spatial distribution of ions, with a 250 μm scale bar.
- グラフ (Graph):** A mass spectrum plot showing intensity versus m/z, with a prominent peak at m/z 799.521.
- MSイメージ一覧 (MS Image List):** A list of MS images, with the selected image highlighted in a green box.
- 解析パラメータ (Analysis Parameters):** A table of analysis settings.

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