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The screenshot displays the IMAGEREVEAL software interface. The 'File' menu is open, and 'Reference Image Management' (参照画像の管理) is highlighted with a yellow box and a mouse cursor. The main window shows a data table with columns for 'No.', '使用', 'タグ', and 'ラベル'. Below the table is a 'グラフ' (Graph) section with a mass spectrum plot titled 'Testicle\_9AA\_PI\_SL\_5x\_1\_AREA01.imdx Whole\_Ave.' showing intensity versus m/z. The plot has several peaks labeled with their m/z values: 721,482, 744,540, 767,492, 795,521, 796,524, 797,524, 811,514, 837,539, and 885,538. To the right of the graph is an 'MSイメージ一覧' (MS Image List) section showing a thumbnail of the reference image and its TIC (Total Ion Chromatogram) plot.

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The screenshot displays the 'Image Registration' (画像レジストレーション) software interface. The window title is '画像レジストレーション' and it includes standard window controls (minimize, maximize, close). The interface is divided into three main sections: '01. 参照画像の選択' (Reference Image Selection), '02. 粗い画像レジストレーション' (Coarse Image Registration), and '03. 細かい画像レジストレーション' (Fine Image Registration). The '02. 粗い画像レジストレーション' section is currently active, indicated by a blue arrow. The central area shows a large, colorful image with a semi-transparent registration box overlaid on it. The box is tilted and has four corner handles for adjustment. To the left of the image is a 'パラメータ調整' (Parameter Adjustment) panel with directional buttons (up, down, left, right) and a central 'MS' button. Below these buttons, the current position and rotation are displayed: '位置 X : 1036.0', '位置 Y : 1071.7', and '回転 : 21.0'. Below the parameter adjustment is a 'ズーム' (Zoom) panel with sliders for X and Y axes. At the bottom left, there is a '調整値の細かさ' (Adjustment Value Fineness) panel with two radio buttons: '細かい' (Fine) and '粗い' (Coarse), with '粗い' selected. To the right of the image is a '参照画像の設定' (Reference Image Settings) panel with sliders for '輝度' (Brightness), 'コントラスト' (Contrast), and '透明度' (Opacity), and a dropdown for '平滑化フィルター' (Smoothing Filter) set to 'なし' (None). Below this is an 'MSイメージ設定' (MS Image Settings) panel with a 'ファイル' (File) dropdown and an 'MSイメージ' (MS Image) dropdown set to 'TIC'. At the bottom of the window, there is a color scale for the 'MSイメージ' and a set of navigation buttons: '戻る' (Back), '次へ' (Next), '保存' (Save), and 'キャンセル' (Cancel).

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The screenshot displays the 'Image Registration' (画像レジストレーション) software interface. The process is divided into three stages: 01. Reference Image Selection, 02. Rough Image Registration, and 03. Fine Image Registration (highlighted in blue). The central workspace shows a satellite image with a red grid overlay, indicating alignment. The left sidebar contains 'Parameter Adjustment' (パラメータ調整) with 'Manual' (手動) selected, a 'Show Grid' (グリッドを表示する) checkbox, a grid size slider, and an 'Automatic Registration Execution' (自動レジストレーション実行) button. The right sidebar shows 'Reference Image Settings' (参照画像の設定) with sliders for 'Brightness' (輝度), 'Contrast' (コントラスト), and 'Transparency' (透明度), and a 'Smoothing Filter' (平滑化フィルター) dropdown set to 'None' (なし). At the bottom, there are navigation buttons: 'Back' (戻る), 'Next' (次へ), 'Save' (保存), and 'Cancel' (キャンセル). A hand cursor is pointing at the 'Save' button.