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The screenshot displays the main interface of the IMAGEREVEAL software. The window title is "差異解析 | No Title\* - IMAGEREVEAL". The interface is organized into several functional areas:

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The screenshot shows the MS Image software interface. At the top, there is a toolbar with several icons and a color scale. Below the toolbar is a large image area displaying a colorful, noisy pattern. To the right of the image area is a panel with text information. At the bottom right, there is a button labeled '情報のコピー'.

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The screenshot displays the software's main interface. On the left is a sidebar with icons for file operations, image settings, and data processing. The top menu bar includes '設定' (Settings), 'オプション' (Options), '拡張機能' (Extensions), and 'ヘルプ' (Help). A context menu is open over the '設定' button, with 'MSイメージ...' highlighted. The main window is divided into several panels: a data table at the top, a mass spectrum graph at the bottom left, and an MS image viewer at the bottom right. The mass spectrum graph shows intensity versus m/z with several peaks labeled. The MS image viewer shows a color-coded image of the sample with a scale bar of 250 μm.

No.	Name
	Testicle_9AA_Pi_SL_5x_1_AREA01.imdx Whole_Ave.

m/z	Intensity
721.4819	
744.5405	
767.4918	
795.5208	
796.5236	
797.5237	
811.5141	
837.5390	
885.5378	

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The screenshot displays the software's main interface. The top menu bar includes '設定' (Settings), which is highlighted with a yellow box and a hand cursor. A dropdown menu is open, showing options like 'アダクトイオン...', '化合物テンプレート...', 'マトリックス...', 'MSイメージ...', 'ROI ...', 'グラフ...', and 'データ行列・ストレージ...'. The 'ROI ...' option is highlighted with a yellow box and a hand cursor. The main window is divided into several panels: a sidebar on the left with various tool icons, a central data table with columns for 'No.', 'ファイル名', 'ROI名', and '属性', a 'グラフ' (Graph) panel showing a mass spectrum plot for 'Testicle\_9AA\_PI\_SL\_5x\_1\_AREA01.imdx Whole\_Ave.' with peaks labeled at m/z values such as 721,482, 747,495, 767,492, 794,524, 795,521, 810,512, 837,539, 857,509, and 885,538, and a 'MSイメージ一覧' (MS Image List) panel at the bottom right showing a heatmap of the sample.

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