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The screenshot displays the IMAGEREVEAL software interface. A central dialog box titled '解析手法の選択' (Analysis Method Selection) is open, showing a list of analysis methods under the heading '新規プロジェクト' (New Project). The methods listed are:

- 差異解析 (Differential Analysis)
- 画像解析 (Image Analysis)
- 定量解析 (Quantitative Analysis) - This option is highlighted with a yellow border and a mouse cursor is pointing at it.
- スクリーニング (Screening)
- スクリーニング (定量) (Screening (Quantitative))

The background interface includes a sidebar with various analysis tools, a main workspace with data tables and graphs, and a right-hand panel with sample information. The top menu bar shows 'ファイル', '表示', '検索', '設定', 'オプション', '拡張機能', and 'ヘルプ'. The title bar indicates '差異解析 | No Title - IMAGEREVEAL'.

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The screenshot displays the '定量解析' (Quantitative Analysis) window in the IMAGEREVEAL software. The interface is divided into several panels:

- Left Panel (Navigation):** A vertical sidebar with icons for various functions. The 'IMDXファイルを追加' (Add IMDX File) icon is highlighted with a yellow box and a mouse cursor.
- ROIリスト (ROI List):** A table with columns: No., 使用 (Use), ファイル名 (File Name), ROI名 (ROI Name), and 属性 (Attribute).
- 検量線 (Calibration Curve):** A graph showing Concentration on the y-axis (0 to 10) and Intensity on the x-axis (0 to 10). A checkbox '原点を通る' (Pass through origin) is present.
- Parameters:** Input fields for '許容幅' (Tolerance) set to 0.2000 Da and 'しきい値' (Threshold) set to 0.00000. A button 'MSイメージ一覧に追加' (Add to MS Image List) is also visible.
- MSイメージ (MS Image):** A large dark area for the MS image, with a '情報のコピー' (Copy Information) button.
- MSイメージ一覧 (MS Image List):** A panel at the bottom right with controls for adding, deleting, and sorting MS images, and buttons for '重ね合わせ' (Overlay) and '四則演算' (Arithmetic).
- ROI平均 (ROI Average):** A small table with columns '計..' (Calculation), '表..' (Table), and 'ファイル名' (File Name).
- 解析パラメータ (Analysis Parameters):** A section for 'ノーマライズ未計算' (Normalization not calculated) with a table for No., Name, and Value.

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The screenshot displays the IMAGEREVEAL software interface with the following components:

- Left Panel:** A vertical menu with options like 'ファイル', '画像設定', 'ROI設定' (highlighted with a yellow box), '前処理設定', '前処理', '検査線', '対象リスト設定', and '検査線計算'.
- ROI List Table:** A table with columns 'No.', '使用', 'ファイル名', 'ROI名', and '属性'. It lists 9 items, all with 'Group A' as the attribute. The table is highlighted with a green border.
- Measurement Line Table:** A table with columns 'No.', 'm/z', and '化合物名', currently empty.
- Concentration vs Intensity Plot:** A graph with 'Concentration' on the y-axis (0-10) and 'Intensity' on the x-axis (0-10).
- MS Image:** A large heatmap image with a color scale and a '250 μm' scale bar. Metadata on the right includes '化合物名/コメント: TIC', 'ファイル名: Slide2-2\_CHCA蒸着\_Sec2\_5x\_10um\_AREA01.imdx', and 'タイプ: TIC'.
- MS Image List:** A grid of smaller MS images, each labeled 'TIC'.
- Mass Spectrometry Plots:** Four stacked plots showing 'Intensity' vs 'm/z' for different ROI areas (AREA01 to AREA04). The x-axis for all plots ranges from 610 to 660 m/z.
- ROI Average Table:** A table with columns '計..', '表..', and 'ファイル名'.
- Analysis Parameters:** A section titled '解析パラメータ' with a 'ノーマライズ未計算' button and an empty table with columns 'No.', 'Name', and 'Value'.

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
The screenshot displays a software interface for quantitative imaging analysis. The main window is titled "定量解析" (Quantitative Analysis) and shows the following components:

- ROIリスト (ROI List):** A table listing 17 ROIs. The "使用" (Use) column has checkmarks for ROIs 9 through 17. The "ROI名" (ROI Name) column shows "All" for ROIs 1-8 and "ROI001" through "ROI008" for ROIs 9-17. The "属性" (Property) column is "Group A" for all.
- 検量線 (Calibration Curve):** A graph showing Concentration vs. Intensity. The y-axis ranges from 0 to 10, and the x-axis ranges from 0 to 10. A checkbox "原点を通る" (Pass through origin) is present.
- MSイメージ (MS Image):** A large heatmap showing the intensity distribution of the sample. A green circle highlights a specific region labeled "ROI004".
- MSイメージ一覧 (MS Image List):** A grid of smaller heatmaps showing the intensity distribution for different ROIs. The first row shows four heatmaps labeled "180406\_Right1-2\_Right1\_1.25x+Matrix\_BW\_AREA01.imdx Whole\_Ave.", "180406\_Right1-2\_Right1\_1.25x+Matrix\_BW\_AREA02.imdx Whole\_Ave.", "180406\_Right1-2\_Right1\_1.25x+Matrix\_BW\_AREA03.imdx Whole\_Ave.", and "180406\_Right1-2\_Right1\_1.25x+Matrix\_BW\_AREA04.imdx Whole\_Ave.". The second row shows four heatmaps labeled "180406\_Right1-2\_Right1\_1.25x+Matrix\_BW\_AREA05.imdx Whole\_Ave.", "180406\_Right1-2\_Right1\_1.25x+Matrix\_BW\_AREA06.imdx Whole\_Ave.", "180406\_Right1-2\_Right1\_1.25x+Matrix\_BW\_AREA07.imdx Whole\_Ave.", and "180406\_Right1-2\_Right1\_1.25x+Matrix\_BW\_AREA08.imdx Whole\_Ave.". Each heatmap is labeled "TIC" (Total Ion Chromatogram).
- 対象リスト設定 (Target List Setting):** A sidebar menu on the left with a yellow box around the "対象リスト設定" (Target List Setting) option.
- 解析パラメータ (Analysis Parameters):** A section at the bottom left with a "ノーマライズ未計算" (Normalization not calculated) button and a table with columns "No.", "Name", and "Value".

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Matrix Clusters  
Lipids  
Lipid Mediators  
Amiodarone

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The screenshot displays the IMAGEREVEAL software interface for setting calibration curves. The main window is divided into several panels:

- ROIリスト (ROI List):** A table listing 17 ROIs with their respective file names, ROI names, and attributes (Group A).
- 検量線 (Calibration Curve):** A graph showing Concentration (Y-axis, 0 to 10) versus Intensity (X-axis, 0 to 10). Below the graph are controls for tolerance (許容幅: 0.2000 Da) and threshold (しきい値: 0.00000).
- MSイメージ (MS Image):** A heatmap showing the distribution of ions across the sample area. A red circle highlights the ROI (ROI004). The image includes a 100 μm scale bar and a color scale legend.
- MSイメージ一覧 (MS Image List):** A grid of smaller MS images for different ROIs, each labeled with 'TIC' (Total Ion Chromatogram).

The left sidebar contains various tool icons, with the '検量線計算' (Calibration Curve Calculation) icon highlighted in a yellow box. The bottom panel shows the '解析パラメータ' (Analysis Parameters) section, including a table for 'ノーマライズ未計算' (Not Normalized) parameters.

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