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The screenshot displays the IMAGEREVEAL software interface. The main window is divided into several panels:

- ROIリスト (ROI List):** A table listing 16 ROIs. The 'データ点数' (Data Points) column is highlighted with a green box. A callout box points to this column with the text: "ここに各データ点数が表示されます" (Data points for each ROI are displayed here).
- データ行列テーブル (Data Matrix Table):** A table showing 17 rows of data with columns for No., 使用 (Use), タグ (Tag), ラベル (Label), m/z, and PCA-水平軸 (PCA Horizontal Axis).
- MSイメージ (MS Image):** A mass spectrum image showing a red and white pattern. 15 ROIs are marked with colored circles and labeled (ROI001 to ROI015). A 250 μm scale bar is visible in the bottom right.
- ROI平均 (ROI Average):** A panel for calculating averages for all ROIs, with a list of files and checkboxes.
- MSイメージ一覧 (MS Image List):** A panel showing a list of MS images and their corresponding spectra. The spectra are labeled with EXP = ((A+B)/C) and m/z ranges.
- 解析パラメータ (Analysis Parameters):** A table showing parameters such as normalization, data matrix analysis method, m/z range, and PCA calculation options.

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The interface is divided into several panels:

- ROIリスト (ROI List):** A table listing 16 ROIs. ROI007 is selected, and a context menu is open over it.
- データ行列テーブル (Data Table):** A table showing data points for the selected ROI007.
- MSイメージ (MS Image):** A mass spectrometry image showing the spatial distribution of ions, with ROI007 highlighted in green.

The context menu for ROI007 includes the following options:

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The data table below shows the following data points:

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