

Data conversion Actual result

※not guaranteed

Actual results from imzML to imdx

Vendor	Platform	Software to export imzML
Shimadzu	iMScope TRIO™	DataConverter
Kratos	MALDI-7090™	MALDI Solution™
Bruker	Solalix™, UltraFlex™	FlexImaging™ 4.9, 5.0(Processed, normalized)
Thermo Fisher	Q Exactive™(Orbitrap™)	raw to imzML Converter 1.1.4.5i Beta *
Waters	DESI-Q-TOF	HD Imaging™ V.1.4
Alexandrov Team at EMBL Heidelberg.		pyimzML 1.2.1 **
Alan M. Race, et.al University of Birmingham		imzMLConverter 1.3 ***

* <https://ms-imaging.org/wp/imzml/software-tools/raw-to-imzml-converter/>

** <http://pyimzml.readthedocs.io/en/latest/#>

*** <https://www.cs.bham.ac.uk/~ibs/imzMLConverter/>

Third- party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether not they are used with trademark symbol "TM" or ®".